

Power	220VAC ±10%
Power Consumption	180 Watts
Current	0.8 A
Medication Capacity	6ml
Particle Size	0.5-5µm
MMAD	2.95 µm
Sound Level	<60 dBA
Average Nebulization Rate	Min 0.2ml/min
Compressor Pressure Range	30 to 50 Psi
Operating Pressure Range	9.5 to 16 Psi
Liter Flow Range	6 ~ 10 lpm
Operating Temperature Range	10 °C to 40 °C
Operating Humidity Range	10 to 95% RH
Storage Temperature Range	-20 °C to 70 °C
Storage Humidity Range	10 to 95 % RH
Dimension (L x W x H)	188 x 106 x 188mm
Weight	3.3Lbs (1.5kg)
Standard Accessories	Nebulizer Cup, Connecting Tube, Mouthpiece, Adult & Pedia Mask, 5pcs Replacement Filters

SPECIFICATIONS



Recommended Replacement Kit

Your Compressor Nebulizer should be used under the supervision of a licensed physician and/or a respiratory therapist. We encourage you to thoroughly read this guidebook to learn about the features of this product. Any use of this product other than its intended use should always be avoided.

This product is developed for the effective treatment of asthma, allergies and other respiratory disorders. By creating a stream of air that travels through the clear air tube to the nebulizer. Once the steam of air enters the nebulizer, it will convert the prescribed medication into an aerosol mist for easy inhalation.


This product is designed to efficiently deliver physician prescribed medication to bronchial lung passages. With proper care and use, this unit will provide you with years of reliable treatment.

INTRODUCTION

Thank you for purchasing TMS Compressor Nebulizer. The compressor nebulizer is a compact medical device designed to efficiently deliver physician prescribed medication to bronchial lung passages. With proper care and use, this unit will provide you with years of reliable treatment.

Your Compressor Nebulizer should be used under the supervision of a licensed physician and/or a respiratory therapist. We encourage you to thoroughly read this guidebook to learn about the features of this product. Any use of this product other than its intended use should always be avoided.

PRODUCT IDENTIFICATION



IMPORTANT SAFEGUARDS

The following basic precautions are needed when using electrical products:

Cautions: Failure to read and observe all precautions could result in personal injury or equipment damage.

1. To avoid electrical shock: Keep the unit away from water - Do not immerse the power cord or the unit in liquid - Do not use while bathing - Do not reach for a unit that has fallen into water - Immediately unplug the unit.

2. Never operate the unit if it has damaged parts (including power cord), if it has been dropped or submerged in water. Promptly send it to a service center for examination and repair.

3. The unit should not be used where flammable gas, oxygen or aerosol spray products are being used.

4. Disconnect the unit from the electrical outlet before cleaning, filling and after use.

5. Do not use other attachments unless recommended by manufacturer.

The compressor motor has a thermal protector which will shut off the unit before the unit is overheated. When the thermal protector shuts the unit off, please:

- Switch off the unit.
- Unplug the unit from the electrical outlet.
- Wait 30 minutes for the motor to cool down before another treatment. Make sure the air openings are not obstructed.

Important:

- When treatment is finished, shut off the unit and unplug it from the electrical outlet.
- Press power switch to begin your prescribed treatment.
- Plug the power cord into an appropriate electrical outlet. Make sure at this stage, the power switch is at "OFF" status.
- Attach mouth piece or mask per your choice to the top section of nebulizer.
- Attach the other end of air tube to the air tube connector located on the front of the compressor.
- Attach one end of air tube to the base of the nebulizer sections fit well.
- Gently twist the top part clockwise to reassemble the nebulizer. Be sure the two physician. Be sure that the cone is put inside the bottom section.
- Fill the bottom section of the nebulizer with the medication prescribed by your physician. Be sure that the cone is put inside the bottom section.
- Gently twist the top part of the nebulizer counter-clockwise to disassemble the nebulizer.
- Open the top cover.
- Place your Compressor Nebulizer on a flat and stable surface. Be sure that you can easily reach the controls when you are seated.
- Clean all necessary parts after each use as instructed in this guidebook.

OPERATING YOUR COMPRESSOR NEBULIZER

- Do not immerse the unit in water. It may damage the unit.
- Disconnect the unit from the electrical outlet before cleaning.
- Clean all necessary parts after each use as instructed in this guidebook.

CLEANING CAUTIONS


- Do not store the unit under direct sunlight, high temperature or humidity.
- Always keep the unit unplugged while not in use.
- Keep the unit out of reach of small children.

STORAGE CAUTIONS

- Connect this product to the appropriate voltage outlet for your model.
- Do not operate this product then be left unattended.
- Never operate if this unit has a damaged cord or plug, if it has been dropped into water, or is not working properly, return the unit to a service center for repair.
- If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
- Always unplug the product immediately after use.
- Never block the air openings of the main unit or place it where air openings may be obstructed.

OPERATING CAUTIONS

- Connect this product to the appropriate voltage outlet for your model.
- Do not operate this product then be left unattended.
- Never operate if this unit has a damaged cord or plug, if it has been dropped into water, or is not working properly, return the unit to a service center for repair.
- If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
- Always unplug the product immediately after use.
- Never block the air openings of the main unit or place it where air openings may be obstructed.



Please read the instructions

MODEL: T200

- Do not use cotton or any other materials. Do not wash or clean the filter. Only use filters supplied by your distributor, and do not operate without a filter.
- Change the filter every 30 days or when the filter turns gray.
- Changing procedure:
 - Remove the filter cover
 - Replace the used filter with a new one.
 - Replace the filter cover.

FILTER CHANGE

- Wipe daily with a damp cloth.
- Do not use any powdered cleaners or soap pads, which may damage the finish.

CLEANING THE COMPRESSOR

- Use one part white vinegar with three parts distilled water. Make sure mixed solution is enough to submerge the nebulizer, mouthpiece and mask.
- Complete Rinsing Steps 1-3.
- Wash nebulizer, mouthpiece and mask in warm water and a mild detergent. Then wash them in hot tap water.
- Submerge these parts in the vinegar and water solution for thirty minutes.
- Complete Rinsing Steps 3-5.

DISINFECTIONS

- Disconnect the air tube, nebulizer, mouthpiece and mask.
- Gently twist the nebulizer to open it.
- Rinse the nebulizer, mouthpiece and mask with water.
- Dry them with clean soft towel or let it air dry.
- Reassemble the nebulizer when completely dry and put these parts in a dry, sealed container.

RINSING (after each treatment)

It is recommended that the nebulizer, mouthpiece and mask are thoroughly cleaned with hot water after each use and cleaned with a mild detergent after the last treatment of the day. If your physician or respiratory therapist specifies a different cleaning procedure, express to them clearly your nebulizer model.

CLEANING

Caution: Before cleaning, make sure the power cord is disconnected with the main outlet and the unit.



Product Warranty Card

NAME:	DATE OF PURCHASE:
ADDRESS:	PURCHASE FROM:
ITEM PURCHASED: TMS NEBULIZER PUMP	REGISTRATION DATE:

Inmed Corporation warrants this product to be free from defects in material or workmanship within the specified warranty period under normal use. If fault is found, please return the equipment to the store where product was purchased. Inmed Corporation will repair or replace any defective part free of charge subject to the terms and conditions stated herein.

For service, the unit is to be returned freight prepaid to:

Inmed Corporation

5 Calle Industria, Bagumbayan,
Quezon City 1110, Philippines
Tel: 02.8571.1888

Please register your unit online at www.inmed.com.ph