




anchor

surgical needles
and other OR related products



Authorized Distributor:

INMED
CORPORATION

Tel.: +63.2.5711888 • Fax: +63.2.5719912
www.inmed.com.ph



contents

SOFT TISSUE BIOPSY DEVICES		4
OPHTHALMIC BLADES		5
SHARPS COUNT SYSTEMS		6
Count-Tainer® 20 (foam)	5120F	
Count-Tainer® 20 (magnetic)	5220M	
Count-Tainer® 40 (foam)	5100F	
Count-Tainer® 40 (magnetic)	5200M	
Count-Tainer® 50 (foam block)	5350	
OTHER O.R. PRODUCTS		7
Trocar Port Suture Needle	TR-1001	
Myringotomy Blade (with handle)	MB-3003	
Reusable stainless sterilizing rack		
Reusable stainless sterilizing box		
NEEDLE PACKAGING		7
CUSTOM NEEDLE KITS		8
SURGICAL NEEDLES	Style No.	Page
Abscission	1851	9
Bunnels	1854	9
Cervix	1835	9
Dura	(see 1833 and 1842)	9
Eye	1830, 1831, 1833 1857	9
Ferguson	1842	9
Fish Hook	1853	10
Fistula	1832	10
French Spring Eye	1861, 1862	10
Hagedorn	1852	10
Hernia	1829, 1847	10
Intestinal	1843, 1850	10
Keith	1827	11
Kidney	1838, 1839, 1840	11
King	1841	11
Lanes, Cleft Palate	1844, 1848	11
Loopuyt's	1846	11
Martin Uterine	1860	12
Mayo	1824	12
Mayo Intestinal	1823, 1828, 1837	12
McGowan	1859	12
Murphy Intestinal	1825	12
Orthopedic	1869, 1870	12
Regular Surgeons	1821 1822, 1834,	13
Tonsil, Davis	1845, 1849	13
Trocar	1826, 1836	13
DENTAL NEEDLES		14
POST MORTEM & VETERINARY NEEDLES		15 and 16
NEEDLE CONVERSION CHART: Sterile Disposable		17
NEEDLE CONVERSION CHART: Non-Sterile Reusable		18
VETERINARY SUTURE CASSETTES		19

Authorized Distributor:



Tel.: +63.2.5711888 • Fax : +63.2.5719912
www.inmed.com.ph



anchor

SOFT TISSUE BIOPSY DEVICES



ONE HAND OPERATION

- innovative, comfort-designed, pistol-grip handle fits any hand
- permits excellent control during soft tissue biopsy procedures
- automatic locking assures secure, consistent, accurate sample removal
- sharp cutting needle of polished stainless steel for smooth penetration
- single-use, disposable, sterile device, with several gauges and lengths to choose from
- economical, competitively priced, American made

ORDERING INFORMATION:

SPECIFY CODE NO.

Code No.	Gauge	Notch	Length	
ABD6-14	14	19mm	6"	15.2cm
ABD45-14	14	19mm	4½"	11.4cm
ABD6-18	18	19mm	6"	15.2cm
ABD45-18	18	19mm	4½"	11.4cm
ABD6-20	20	19mm	6"	15.2cm
ABD3-20	20	19mm	3"	7.6cm
ABD15-14*	14	19mm	15"	38.1cm
ABD15-18*	18	19mm	15"	38.1cm
ABD45-14S	14	9.5mm	4½"	11.4cm
ABD3-20S	20	9.5mm	3"	7.6cm

EACH DEVICE PACKAGED STERILE/DISPOSABLE
AVAILABLE IN BOXES OF 5 UNITS OR CASES OF 20 UNITS

*Available in single unit cartons only

Authorized Distributor:

INMED
CORPORATION

Tel.: +63.2.5711888 • Fax: +63.2.5719912
www.inmed.com.ph



AB101



AB102




AB103



AB104




AB200 Series



AB300 Series



AB401



AB500 Series

PRODUCT CODE	DESCRIPTION	PACKAGING
--------------	-------------	-----------

Stab Incision Eye Knives

AB101	Blade Length 5.7mm, Blade Angle 15°	6 per box
AB102	Blade Length 5.7mm, Blade Angle 22.5°	6 per box
AB103	Blade Length 5.7mm, Blade Angle 30°	6 per box
AB104	Blade Length 5.7mm, Blade Angle 45°	6 per box

Phaco/Slit Knives

AB201S	Blade Width 2.5mm, Straight	6 per box
AB201A	Blade Width 2.5mm, Angled	6 per box
AB202S	Blade Width 2.8mm, Straight	6 per box
AB202A	Blade Width 2.8mm, Angled	6 per box
AB203S	Blade Width 3.0mm, Straight	6 per box
AB203A	Blade Width 3.0mm, Angled	6 per box
AB204S	Blade Width 3.2mm, Straight	6 per box
AB204A	Blade Width 3.2mm, Angled	6 per box
AB205S	Blade Width 3.5mm, Straight	6 per box
AB205A	Blade Width 3.5mm, Angled	6 per box

IOL Round Tip Implant Knives


AB301S	Blade Width 4.0mm, Straight	6 per box
AB301A	Blade Width 4.0mm, Angled	6 per box
AB302S	Blade Width 5.2mm, Straight	6 per box
AB302A	Blade Width 5.2mm, Angled	6 per box

MVR Corneal/Scleral Blade

AB401	20 Gauge	6 per box
-------	----------	-----------

Crescent Knives

AB501S	Straight	6 per box
AB501A	Angled	6 per box



IMPROVED COUNT-TAINER®, INJECTION MOLDED IN TOUGH IMPACT PLASTIC

Count-Tainer® 20

Developed for minor procedures, as a compact version of the widely accepted Count-Tainer® 40, complete with friction-fit slotted foam.*

Friction-fit slotted foam provides for a more rapid technique—saves time by just popping the needle in place in an ordered format—no need to reclamp needles.

Non-flaking foam strips minimize particulate matter.

Large, legible numerals for ease in counting.

Backing tape secures to any surface—a “no hands on product” procedure.

Rigid plastic construction with positive lock eliminates needle punctures and escaping sharps.

Adhesive well for larger needles, blades, and miscellaneous sharps.

Available in both a slotted foam version (Code number 5120F) and a magnetic version (Code number 5220M). Both versions are pouched sterile, disposable—20 units to a dispenser, 80 units to a case (4 dispensers).

Count-Tainer® 40


This widely accepted sharps count system is intended for use in major procedures. Count-Tainer® 40 has the identical features and benefits described above in the Count-Tainer® 20.

Available in both a slotted foam version (Code number 5100F), and a magnetic version (Code number 5200M). Both versions are pouched sterile, disposable—10 units to a dispenser, 40 units to a case (4 dispensers).

Count-Tainer® 50

Large foam block design for major procedures. Foam block ensures secure placement, with room for 50-100 needles. Adhesive well for larger needles, blades and miscellaneous sharps. Pouched sterile, disposable—64 units to a case. CODE NUMBER 5350.


* Patent No. 4243140



TROCAR PORT SUTURE NEEDLE

NEW, from ANCHOR, a TROCAR PORT SUTURE NEEDLE for use in laparoscopic surgery. A Simple, Effective and Economical way to reduce Trocar Port Hernias. It has a cutting edge point for penetrating fascia and an ergonomic handle for superior control.

Each box contains one needle in a sterile pouch.
Product Code #TR-1001



MYRINGOTOMY "DOUBLE EDGE" BLADE

Sterile Disposable MYRINGOTOMY BLADE, with a 30 degree bend, for use in E.N.T. surgery. The blade is manufactured from quality stainless steel, 20 gauge diameter, which is ground, honed and 100% inspected to ensure that the edge meets our high performance standards.

Each box contains 20 sterile single packed blades.
Sterile with handle (length: 7") Product Code #MB-3003



NEEDLE STERILIZING RACK

Reusable stainless steel rack with spring holder, firmly aligns up to 24 surgical eyed needles during sterilization.

Available in three models: the Narrow Spring Holder, #1700, for smaller needles; the Straight Spring Holder, #1800, for large gauge needles; and the Tapered Spring Holder, #1900, for a variety of small and large needles.

NEEDLE STERILIZING BOX

Reusable two-piece stainless steel container designed to accommodate the Needle Sterilizing Rack and other miscellaneous stainless products during sterilization.

Sterilizing Box 4" x 5" x 1/2"	No. 1798
Narrow Spring Holder	No. 1700
Straight Spring Holder	No. 1800
Tapered Spring Holder	No. 1900

Combinations of Boxes, Holders and Needles:

Box with No. 1700 Holder	No. 1798-17
Box with No. 1800 Holder	No. 1798-18
Box with No. 1900 Holder	No. 1798-19

Note: The above spring holders are available with 2 dozen needles — simply add on an A after code number. (Example: 1700A)



NEEDLE PACKAGING

Anchor *Disposable* Needles are packed in units of two. Each pair of needles is double-pouched in "see-through," peel-apart, sterile Tyvek®.

They are available in small cartons of 40 needles (20 pouches), and a dispenser-pack gross of 144 needles (72 pouches). Order by the needle number and add the suffix "DC" for disposable cartons, example: 1824-4DC. Add the suffix "DG" for disposable gross boxes, example: 1824-4DG.

Anchor *Reusable* Needles are non-sterile and are packed in units of six (minimum order: 2 packages, 1 dozen). Each unit is contained in the widely known Anchor purple glassine envelope and sleeve that clearly identifies the quantity, style and size. Order by the needle number and add the suffix "R" for reusable, example: 1824-4R.



PROPER NEEDLE USE

Improper surgical technique in the use of a suturing needle may result in needle breakage. Small needle holders should be used on fine gauge needles and heavy holders on heavy gauge needles to avoid needle breakage because of undue stress. The weakest part of a needle is the eye and the point. Anchor needles are manufactured with a flat area on the shank to facilitate proper clamping with holder and to prevent turning in holder. The needle holder should be applied on the needle flat part about 1/4 of length of needle from the eyed end. **DO NOT CLAMP HOLDER DIRECTLY ON THE EYE, OTHERWISE BREAKAGE MAY OCCUR.** When withdrawing needle, clamp needle on shank about 1/4 length from the point. **DO NOT GRASP NEEDLE POINT, OTHERWISE BREAKAGE MAY OCCUR.** "Examine needle carefully after each use. If needle becomes bent or damaged during use, withdraw and discard needle. **DO NOT RE-USE BENT OR DAMAGED NEEDLE, OTHERWISE BREAKAGE MAY OCCUR.**"



ANCHOR-PAK...the time saver; the money saver.

The ultimate in eyed-needle organization. Take the time now to set up *only once*, your needle requirements for any given surgical procedure. Send the set-up to ANCHOR and we will duplicate the set-up in a kit format. Each individual needle will be aligned in a unique friction-fit slot on a foam strip. Since there is no needle penetration of the foam, the risk of particulate matter is minimized. Simply lift out as needed.

ANCHOR will also custom imprint the kit by hospital name, procedure, and number of needles. Each kit is double-pouched in Tyvek®, sterilized, and shipped—organized and ready for use.

It's as good as it sounds. So, take the time now to do your set-up, and send it to Anchor.

We have a feeling it may be *your last set-up.*

Ask about our Sample Kit Program.

SURGICAL NEEDLES



Needles are available both sterile (disposable) and non-sterile (reusable) unless otherwise indicated.

Anchor Disposables: Order by needle number and add the suffix "DC" for disposable cartons of 40 needles (20 pouches), example 1824-4DC. Add the suffix "DG" for disposable gross boxes of 144 needles (72 pouches), example: 1824-4DG.

Anchor Reusables: Order by needle number and add the suffix "R" for reusable, example: 1824-4R.

All needles are shown full scale in this catalog. For your convenience a cross-reference guide is included to show styles and sizes equivalent to other brands of surgical needles.



Authorized Distributor:



Tel.: +63.2.5711888 • Fax : +63.2.5719912
www.inmed.com.ph



ABSCISSION $\frac{3}{8}$ circle

1851




BUNNELS straight

1854




CERVIX $\frac{3}{8}$ circle


1835



EYE



$\frac{3}{8}$ circle
1830



$\frac{1}{2}$ circle
1831



$\frac{1}{2}$ circle
1833



(1833-3, 1833-4 also used for dura)



Jameson Strabismus
1857



(non-sterile only)



FERGUSON $\frac{1}{2}$ circle

1842


(1842-20, 1842-22 also used for dura)







FISH HOOK

1853-3





HERNIA



(1824-6, 1849 also used for hernia)


FISTULA 1/2 circle
(heavy gauge)

1832





INTESTINAL straight
(fine)

1843




FRENCH SPRING EYE 1/2 circle

1861




1862





INTESTINAL 1/2 circle
(fine)

1850



HAGEDORN straight
(flat body)


1852
(non-sterile only)



All needles are shown full scale in this catalog.


KEITH straight

1827





KING 3/8 circle



1841



KIDNEY




LANES CLEFT PALATE




LOOPUYT'S

1846
(non sterile only)







MARTIN UTERINE ½ circle 1860
(heavy gauge) (also used for orthopedics)




MAYO INTESTINAL straight 1828




MAYO ½ circle 1824




McGOWAN ½ circle 1859
(heavy gauge)




MAYO INTESTINAL ½ circle 1823



MURPHY INTESTINAL ½ circle 1825



MAYO INTESTINAL ¾ circle 1837



ORTHOPEDIC






All needles are shown full scale in this catalog.

REGULAR SURGEON'S 1/2 curved (non-sterile only)

1821




CUTTING
EDGE



REGULAR SURGEON'S 3/8 circle

1822



CUTTING
EDGE

TONSIL, DAVIS



3/8 circle
1845




TAPER
POINT



1/2 circle
1849

TROCAR 1/2 circle


1826



TROCAR
POINT

REGULAR SURGEON'S 1/2 circle (light gauge)


1834



CUTTING
EDGE

TROCAR 3/8 circle


1836



TROCAR
POINT



½ curved
(non-sterile only)

1821





¾ circle

1822




½ circle

1832




½ circle

1834




½ circle

1844



¾ circle

1848





Post Mortem and Veterinary needles are available non-sterile only

All needles are shown full scale


DOUBLE CURVED POST MORTEM

1819




HALF CURVED VETERINARY

1820




CRUCIATE REPAIR "J" CURVED

1816




VETERINARY straight

1817





VETERINARY 3/8 circle

1818




FISTULA 1/2 circle

1832





REGULAR SURGEON'S 3/8 circle

1822





VETERINARY half curved

1821



LOOPUYT'S

1846



ANCHOR STERILE DISPOSABLE CONVERSION TABLE



This table is intended for use as a guide in converting from other manufacturer's needle code numbers to ANCHOR equivalencies.

There could be instances when an exact conversion on a specific needle is not possible. In these instances, ANCHOR will provide comparable needle styles and let you decide your preference. ANCHOR also has the capability of duplicating virtually all special needles.

For the most accurate conversion, it is recommended that product samples be forwarded.

ANCHOR	RICHARD-ALLAN	HMS - TORRINGTON	ANCHOR	RICHARD-ALLAN	HMS - TORRINGTON	ANCHOR	RICHARD-ALLAN	HMS - TORRINGTON
1821-1	X	X	1831-1	2123-1	709-1	1843-1½	2147-4	721-1½
1821-2	X	X	1831-2	2123-2	709-2	1843-1¾	2147-3	721-1¾
1821-4	X	X	1831-3	2123-3	709-3	1843-2¼	2147-1	721-2¼
1821-6	X	X	1831-4	2123-4	709-4	1843-2¾	X	721-2¾
1821-8	X	X	1832-00	X	725-00	1844-1	2154-1;2157-1	729-1
1821-10	X	702-10	1832-0	2128-1	725-0	1844-2	2154-7/8;2157-2	729-2
1821-12	X	702-12	1832-1	2128-2	725-1	1844-3	2154-6/8;2157-3	729-3
1821-14	X	702-14	1832-2	2128-3	725-2	1844-4	2154-5/8;2157-4	729-4
1821-16	X	702-16	1832-3	2128-4	725-3	1844-5	2154-4/8,3/8	729-5
1821-18	X	702-18	1832-4	2432-4	725-4	1845	2195-0	732
1821-20	X	702-20	1832-5	2128-5;2432-5	725-5	1846-1	X	X
1821-22	X	X	1832-6	2128-6;2432-6	725-6	1846-2	X	X
1822-1	2090-1,2	700-1	1832-7	2432-7	725-7	1846-3	X	X
1822-2	2090-3	700-2	1832-8	2432-8	725-8	1846-4	X	X
1822-4	2090-4,5	700-4	1832-9	2432-9	725-9	1846-5	X	X
1822-6	2090-6,7	700-6	1833-1	2125-1;2500-0	716-1	1847	X	733
1822-8	2090-8,9,10	700-8	1833-2	2092-1;2125-2	716-2	1848-1	2153-1;2158-1	730-1
1822-10	2090-11,12	700-10	1833-3	2092-2;2125-3	716-3	1848-2	2153-7/8;6/8;2158-2	730-2
1822-12	2090-13,14	700-12	1833-4	2092-3;2125-4	716-4	1848-3	2153-5/8;2158-4	730-3
1822-14	2090-15,16	700-14	1834-00	2091-2,2091-3;2127-1	X	1848-4	2153-4/8	730-4
1822-16	2090-17,18,19,20	700-16	1834-0	2091-5	701-0	1848-5	2153-3/8	730-5
1822-18	2090-22,24	700-18	1834-1	2091-6,7;2127-2	701-1	1849	X	731
1822-20	2090-26,28,30	700-20	1834-2	2091-8;2127-3	701-2	1850-1	2150-5	703-1
1823-1	2151-1	705-1	1834-3	2091-9	701-3	1850-2	2150-4	703-2
1823-2	2151-2	705-2	1834-4	2091-10,11,12;2127-4,5	701-4	1850-3	X	703-3
1823-3	2151-3	705-3	1834-5	2091-13,14	701-5	1850-4	2150-3	703-4
1823-4	2151-4	705-4	1834-6	2091-15,16;2127-6	701-6	1850-5	2150-2	703-5
1823-5	2151-5	705-5	1834-7	2091-17,18;2127-7	701-7	1850-6	2150-1	703-6
1823-6	2151-6	705-6	1834-8	2091-19,20,22;2127-8	701-8	1851-1	X	737-1
1824-1	2167-0	707-1	1834-9	2091-26,28,30	701-9	1851-2	X	737-2
1824-2	2167-1;2222-1	707-2	1835-1	X	717-1	1851-3	X	737-3
1824-3	2167-2,3;2222-2	707-3	1835-2	X	717-2	1851-4	X	737-4
1824-4	2167-3,5;2222-3	707-4	1835-3	2168-1	717-3	1852-1½	X	X
1824-5	2167-4;2222-4	707-5	1835-4	2168-2,3	717-4	1852-2	X	719-2
1824-6	2167-5;2222-5	707-6	1835-5	2168-4	717-5	1852-2½	X	719-2½
1825-1	X	718-1	1835-6	X	717-6	1852-3	X	719-3
1825-2	2159-1;2181-1	718-2	1836-1	2169-1,2	710-1	1853	X	X
1825-3	2140-1,2;2159-2,3;2181-2	718-3	1836-2	2169-3	710-2	1854-2½	2200-2½	X
1825-4	2140-4;2181-3,4	718-4	1836-3	2169-4,5	710-3	1854-2¾	2200-2¾	X
1826-0	X	X	1836-4	2169-6	710-4	1854-2	2200-2	X
1826-1	2170-1	726-1	1837-¾	X	727-¾	1857	X	X
1826-2	2170-2;2426-2	726-2	1837-1	X	727-1	1859-1	X	X
1826-3	X	726-3	1837-1¼	X	727-1¼	1859-2	X	X
1826-4	2170-3;2426-4;2701-0	726-4	1837-1¾	X	727-1¾	1859-3	X	X
1826-5	2170-4,5;2426-5	726-5	1838	X	715	1859-4	X	X
1826-6	2170-6;2426-6	726-6	1839	2050-0	711	1859-5	X	X
1827-1½	2134-10	708-1½	1840	X	712	1859-6	2167-6	X
1827-1¾	2134-9	708-1¾	1841-2½	X	713-2½	1860-1	2164-1	714-0
1827-2	2134-8	708-2	1841-2¾	X	X	1860-2	2164-2	714-2
1827-2½	2134-6,7	708-2½	1841-3	X	X	1860-3	2135-1;2164-3	714-3
1827-2¾	2134-4,5	708-2¾	1841-3¼	2145-3	X	1860-4	2135-2;2164-4	714-4
1827-3½	2134-2,3	708-3½	1841-3½	2145-3½;2441-3,5	713-3½	1860-5	2135-3;2164-5	714-5
1827-4	2134-1	X	1841-3¾	2145-4,3¾	X	1860-6	2135-4	714-6
1828-1½	X	704-1½	1842-00	2095-1,2	720-0	1860-7	2164-6,7	714-7
1828-2	2146-3	704-2	1842-1	2095-3,4	720-2	1860-8	X	714-8
1828-2½	X	704-2½	1842-2	2095-5,6,7	720-4	1861-1	2129-1;2002-2;	734S-1
1828-3	2146-1,2	X	1842-4	2095-8	720-6	1861-2	2129-2;2002-3,4	734S-2
1829	X	728	1842-6	2095-9,10	720-8	1861-3	2129-3,4	734S-3
1830-1	2124-3,4	706-1	1842-8	2095-11,12	720-10			
1830-2	2124-5	706-2	1842-10	2095-13,14,15	720-12			
1830-3	2124-6	706-3	1842-12	2095-16,17,18	720-14			
1830-4	X	706-4	1842-14	2442-16;2706-16	720-16			
			1842-16	2442-18;2707-18	720-18			
			1842-18	2442-20	720-20			
			1842-20	2442-22	720-22			

ANCHOR NON-STERILE REUSABLE CONVERSION TABLE



This table is intended for use as a guide in converting from other manufacturer's needle code numbers to ANCHOR equivalencies.

There could be instances when an exact conversion on a specific needle is not possible. In these instances, ANCHOR will provide comparable needle styles and let you decide your preference. ANCHOR also has the capability of duplicating virtually all special needles.

For the most accurate conversion, it is recommended that product samples be forwarded.

ANCHOR	CINCIN- NATI	MILTEX	S. JACK- SON	JORG- ENSEN	ANCHOR	CINCIN- NATI	MILTEX	S. JACK- SON	JORG- ENSEN	ANCHOR	CINCIN- NATI	MILTEX	S. JACK- SON	JORG- ENSEN
1816-2	X	X	X	X	1827-1 1/2	827-1 1/2	431-1 1/2	X	X	1842-00	X	X	E-50	X
1816-3	X	X	X	X	1827-1 3/4	827-1 3/4	431-1 3/4	E-54	X	1842-2	842-2	126-2	E-57	X
1816-4	X	X	X	X	1827-2	827-2	431-2	X	X	1842-4	842-4	126-4	E-56	X
1817-1	X	X	X	X	1827-2 1/2	827-2 1/2	431-2 1/2	E-14	X	1842-6	842-6	126-6	X	X
1817-2	817-2	X	X	X	1827-2 7/8	827-2 7/8	431-2 7/8	X	X	1842-8	842-8	126-8	E-11,24	X
1817-3	817-3	X	X	X	1827-3 1/2	827-3 1/2	431-3 1/2	E-55	X	1842-10	842-10	126-10	X	6037-10
1817-4	817-4	X	X	X	1827-4	827-4	431-4	X	X	1842-12	842-12	126-12	X	6037-11
1818-1	816-5	1281-1	X	X	1828-1 1/2	828-1 1/2	152-1 1/2	X	X	1842-14	842-14	126-14	E-3,10	6037-12
1818-2	816-4	1281-2,3	X	X	1828-2	828-2	152-2	X	X	1842-16	842-16	126-16	E-12	6037-14
1818-3	816-3	1281-4,5	X	X	1828-2 1/2	828-2 1/2	152-2 1/2	X	X	1842-18	842-18	126-18	X	6037-15
1818-4	X	1281-6	X	X	1828-3	828-3	152-3	X	X	1842-20	842-20	126-20	X	6037-16
1819-1	819-1	1295-2	E-69	X	1829	829	200	X	X	1842-22	842-22	126-22	X	X
1819-2	819-2	1295-4	E-70	X	1830-1	830-1	172-1	X	X	1843-1 1/2	843-1 1/2	122-1 1/2	X	X
1819-3	819-3	1295-5	E-71	X	1830-2	830-2	172-2	X	X	1843-1 3/4	843-1 3/4	122-1 3/4	X	X
1819-4	819-4	1295-7	E-72	X	1830-3	830-3	172-3	X	X	1843-2 1/4	843-2 1/4	122-2 1/4	X	X
1820-0	820-0	1291-1	E-66	X	1830-4	830-4	172-4	X	X	1843-2 3/4	843-2 3/4	122-2 3/4	X	X
1820-1	820-1	1291-3	E-67	X	1831-1	831-1	175-1	E-61	X	1844-1	844-1	156-1	X	7031-13
1820-2	820-2	1291-5	X	X	1831-2	831-2	175-2	X	X	1844-2	844-2	156-2	E-36	X
1820-3	820-3	1291-6	X	X	1831-3	831-3	175-3	X	X	1844-3	844-3	156-3	E-29	X
1820-4	X	X	X	X	1831-4	831-4	175-4	X	X	1844-4	844-4	156-4	X	X
1820-5	X	X	X	X	1832-00	832-00	450-2/0	E-48	X	1844-5	844-5	156-5	X	X
1821-1	821-1	141-1	E-68	X	1832-0	832-0	450-1/0	X	X	1845	845	440	X	X
1821-2	821-2	141-2	X	X	1832-1	832-1	450-1	X	X	1846-1	846-1	1274-1	E-63	X
1821-4	821-4	141-3	X	X	1832-2	832-2	450-2	X	X	1846-2	846-2	1274-2	X	X
1821-6	821-6	141-4	X	X	1832-3	832-3	450-3	E-20	X	1846-3	846-3	1274-3	E-64	X
1821-8	821-8	141-5	X	X	1832-4	832-4	450-4	E-19	X	1846-4	846-4	1274-4	X	X
1821-10	821-10	141-6	X	X	1832-5	832-5	450-5	X	X	1846-5	846-5	1274-5	E-65	X
1821-12	821-12	141-7	X	X	1832-6	832-6	450-6	E-39	X	1847	847	201	X	X
1821-14	821-14	141-8	X	X	1832-7	832-7	450-7	X	X	1848-1	848-1	155-1	E-62	X
1821-16	821-16	141-9	X	X	1832-8	832-8	450-8	X	X	1848-2	848-2	155-2	E-15	X
1821-18	821-18	141-10	X	X	1832-9	832-9	450-9	X	X	1848-3	848-3	155-3	X	X
1821-20	821-20	141-11	X	X	1833-1	833-1	176-1	E-33	X	1848-4	848-4	155-4	X	X
1821-22	821-22	141-12	X	X	1833-2	833-2	176-2	X	X	1848-5	848-5	155-5	X	X
1822-1	922-1	140-1	X	X	1833-3	833-3	176-3	E-53,58	X	1849	849	441	X	X
1822-2	922-2	140-2,3	X	7011-2	1833-4	833-4	176-4	X	X	1850-1	850-1	135-1	X	X
1822-4	922-4	140-4,5	X	7011-4	1834-00	934-00	192-1/0,1,2	E-49	7031-2	1850-2	850-2	135-2	X	X
1822-6	922-6	140-6,7	E-38	7011-6	1834-0	934-0	192-3,4,5	X	7031-4	1850-3	850-3	135-3	E-6	X
1822-8	922-8	140-8,9	E-37	X	1834-1	934-1	192-6	X	7031-6	1850-4	850-4	135-4	E-7	X
1822-10	922-10	140-10,11	E-34	7011-8	1834-2	934-2	192-7,8	E-46	X	1850-5	850-5	135-5	E-2,5	X
1822-12	922-12	140-12,13	X	7011-10	1834-3	934-3	192-9,10	X	7031-8	1850-6	850-6	135-6	E-1	X
1822-14	922-14	140-14,15	E-26	7011-12	1834-4	934-4	192-11	E-43	7031-9	1851-1	851-1	178-1	X	X
1822-16	922-16	140-16,17	E-16	7011-13	1834-5	934-5	192-12,13	E-31	7031-10	1851-2	851-2	178-2	X	X
1822-18	922-18	140-18,19	E-41	7011-14	1834-6	934-6	192-14,15	E-30	7031-11	1851-3	851-3	178-3	X	X
1822-20	922-20	140-20	X	7011-16	1834-7	934-7	192-16,17	E-17	7031-12	1851-4	851-4	178-4	X	X
1823-1	823-1	154-1	X	X	1834-8	934-8	192-18,19	E-18	7031-14	1852-1 1/2	X	1095-1 1/2	X	X
1823-2	823-2	154-2	X	X	1834-9	934-9	192-20	E-42	7031-16	1852-2	852-2	1095-2	X	X
1823-3	823-3	154-3	X	X	1835-1	835-1	210-1	X	X	1852-2 1/2	X	1095-2 1/2	X	X
1823-4	823-4	154-4	E-8	X	1835-2	835-2	210-2	X	X	1852-3	X	1095-3	X	X
1823-5	823-5	154-5	X	X	1835-3	835-3	210-3	E-21	X	1853	X	X	X	X
1823-6	823-6	154-6	X	X	1835-4	835-4	210-4	X	X	1854-2	X	X	X	X
1824-1	824-1	112-1	E-60	6037-6	1835-5	835-5	210-5	X	X	1854-2 1/4	X	X	X	X
1824-2	824-2	112-2	X	6037-8	1835-6	835-6	210-6	X	X	1854-2 1/2	854-2 1/2	X	X	X
1824-3	824-3	112-3	E-45	X	1836-1	836-1	114-1	X	X	1857	857	177	X	X
1824-4	824-4	112-4	X	X	1836-2	836-2	114-2	X	X	1859-1	859-1	120-1	X	X
1824-5	824-5	112-5	E-44	X	1836-3	836-3	114-3	E-27	X	1859-2	859-2	120-2	E-47	X
1824-6	824-6	112-6	X	X	1836-4	836-4	114-4	X	X	1859-3	859-3	120-3	X	X
1825-1	825-1	179-1	X	X	1837-3/4	837-3/4	153-3/4	X	X	1859-4	859-4	120-4	X	X
1825-2	825-2	179-2	E-23	X	1837-1	837-1	153-1	X	X	1859-5	X	X	X	X
1825-3	825-3	179-3	E-22	X	1837-1 1/4	837-1 1/4	153-1 1/4	E-25	6017-11	1859-6	X	X	X	X
1825-4	825-4	179-4	E-4,13	X	1837-1 3/4	837-1 3/4	153-1 3/4	X	X	1860-1	860-1	117-1	X	X
1826-0	826-0	119-1/0	X	X	1838	X	435	X	X	1860-2	860-2	117-2	X	X
1826-1	826-1	119-1	X	X	1839	839	436	X	X	1860-3	860-3	117-3	E-59	X
1826-2	826-2	119-2	X	X	1840	840	434	X	X	1860-4	860-4	117-4	X	X
1826-3	826-3	119-3	X	X	1841-2 1/2	841-2 1/2	195-2 1/2	E-51	X	1860-5	860-5	117-5	X	X
1826-4	826-4	119-4	X	X	1841-2 3/4	841-2 3/4	195-2 3/4	X	X	1860-6	860-6	117-6	X	X
1826-5	826-5	119-5	X	X	1841-3	841-3	195-3	E-28	X	1860-7	860-7	117-7	E-40	X
1826-6	826-6	119.6	X	X	1841-3 1/4	841-3 1/4	195-3 1/4	X	X	1860-8	860-8	117-8	X	X
					1841-3 1/2	841-3 1/2	195-3 1/2	X	X	1861-1	871-1	116-1	X	X
					1841-3 3/4	841-3 3/4	195-3 3/4	E-52	X	1861-2	871-2	116-2	X	X
										1861-3	871-3	116-3	E-9	X



- Stronger brake prevents overspin and suture tangles.
- Cassette fits any cassette rack.
- Safety cap prevents cap loss and contamination.
- Anchor Chromic Gut is superior in quality.
- Polymid™ is a new name for the same, superb, coated white nylon.
- Manufactured in the U.S.A.

Veterinary Sterile Suture Cassettes

Polymid™

Order Code	Size	Length
SE-C-40	4/0	100 M
SE-C-30	3/0	100 M
SE-C-20	2/0	100 M
SE-C-10	1/0	75 M
SE-C-01	1	75 M
SE-C-02	2	50 M
SE-C-03	3	50 M
SE-C-04	4	25 M
SE-C-08	8	25 M

Chromic Gut

Order Code	Size	Length
CG-C-40	4/0	50 M
CG-C-30	3/0	50 M
CG-C-20	2/0	50 M
CG-C-10	1/0	50 M
CG-C-01	1	50 M
CG-C-02	2	50 M
CG-C-03	3	50 M
CG-C-04	4	25 M

Authorized Distributor:

INMED[®]
CORPORATION

5 Calle Industria, Bagumbayan,
Quezon City 1110, Philippines

Tel.: +63.2.5711888

Fax : +63.2.5719912

www.inmed.com.ph