

Product Information



simplex[®]

Plain Uncuffed

With High Volume
Low-Pressure Cuff



Endotracheal Tube

Manufactured from the highest quality implantation tested PVC, sterilised and packaged separately for your utmost convenience. Our Endotracheal Tubes are designed with user practicality and patient comfort in mind, to the strictest internationally recognised standards.



Pilot Balloon

Marked, Responsive. Identification and traceability at a glance with prominently displayed tube size and batch number on balloon. Filled balloon provides easy visual pressure confirmation.



Inflation Tube

Smoothly connected, kink resistant. Patient comfort is maximized with smooth connection. Patient safety assured during inflation and deflation through kink resistance.



Inflation Guide

Double ring markings. Clearly marked lines act as a aid to intubation through laryngoscopy.



Cuff

Smooth, Finely welded. Thin, soft cuff ensures effective seal. Heat welded to provide easy cuff to tube transition. Assured atraumatic intubation.



Tubing

Kink-Resistant, Thermosensitive. Softens at body temperature, conforming to contours of respiratory tract, still remaining sturdy enough to ensure patient safety by resisting kinks.



Radiopaque

Splint Resistant. Maintains tube integrity when required length is cut and connector re-inserted. Vital indicator in verification of intubated position as the high quality tubing has clear radiopacity.



Murphy Eye

Smooth. Atraumatic intubation and extubation assured through the polished smoothness of the Murphy eye, an added opening to the tip.



Tip

Bull-Nose. Smooth, short tip minimizes contact with tracheal mucosa, increases the distance from the carina and ensures traumatic intubation.

Simplex offers complete sizes of oral/nasal Murphy cuffed and uncuffed endotracheal tubes which combines quality, dependability and performance with features that set them apart.

With Cuff

Size	Code
5.0	9-IC-SETC50
5.5	9-IC-SETC55
6.0	9-IC-SETC60
6.5	9-IC-SETC65
7.0	9-IC-SETC70
7.5	9-IC-SETC75
8.0	9-IC-SETC80
8.5	9-IC-SETC85

Without Cuff

Size	Code
2.0	9-IC-SETU20
2.5	9-IC-SETU25
3.0	9-IC-SETU30
3.5	9-IC-SETU35
4.0	9-IC-SETU40
4.5	9-IC-SETU45



ISO 9001
2008 SERIES
CERTIFIED

Inmed Corporation
5 Calle Industria, Bagumbayan
Quezon City 1110, Philippines
Tel: +63.2.5711888
Fax: +63.2.5719912
Web: www.inmed.com.ph

The first ISO 9001 certified medical products trading company.

INMED[®]
CORPORATION