

20TH
ANNIVERSARY
CATALOG
2023 EDITION



MEDICAL EQUIPMENTS
OTHER HEALTH CARE PRODUCTS



The first ISO 9001 certified medical products trading company.

INMED
CORPORATION

We help you grow your business.



Who we are

Inmed Corporation is the country's leading and trusted master distributor of quality local and imported medical products and supplies. With the primary focus of providing solid marketing support, Inmed Corporation has established a dependable structure offering significant value to its resellers and dealers in the Philippines and neighboring Asia Pacific Countries.

Committed to help both local and multinational medical partner companies in the Philippines to achieve their target goals. Inmed Corporation continues to work hand-in-hand with them to extend and further expand market reach. The company utilizes its financial strength to purchase high volume supplies and offer its resellers savings at the best attainable cost.

Managed by a team of top professionals in the industry with combined experience of over 60 years in medical products, sales, marketing and distribution, Inmed Corporation boasts of seasoned round-table negotiations, financial packaging and the best business practices benchmarked on international standards providing economical and fundamentally sound solutions to your business needs.

Our Mission

Our mission is to provide quality local and imported medical products at competitive prices and offer professional customer-oriented service.

Our Vision

Our vision is to be the Philippines' leading master distributor of quality medical products – being the preferred supplier of all resellers in the industry and recognized as a strong partner in sustainable healthcare in the nation.

ISO Certified

INMED Corporation is ISO 9000:2015 certified by TUV, compliant with the International Organization for Standardization based in Geneva, Switzerland. It is our commitment to provide our customers with only high quality products, unparalleled service excellence, and reliable timely delivery of all customer orders. Customer satisfaction is our guarantee.



TABLE OF CONTENTS

2 Crank Manual Hospital Bed	13	Nebulizer Pump	11
2-Wheels Rollator	31	Nebulizer with children inspired designs	9
3 Crank Manual Hospital Bed with IV Stand	13	Ordinary Castor	17
3 Function Electric Hospital Bed with IV Stand	14	Overbed Table	16
4-Wheels Rollator	31	Oxygen Concentrator (10L)	12
5 Function Electric Hospital Bed with IV Stand	14	Patient Lifter	18
Aluminum Wheelchair	34	Pedal Exerciser	28
Aluminum Ramp for Wheelchair	36	Pediatric Wheelchair	32
Aluminum Wheelchair (16" width)	34	Physician Weighing Scale	25
Aluminum Wheelchair (18" width)	34	Pocket Mesh Nebulizer	10
Aneroid Sphygmomanometer	22	Portable Defibrillator	7
Anti-Bedsore Mattress System	15	Portable Oxygen Concentrator	11
Baby Weighing Scale	25	Portable Suction Machine	7
Bariatric Wheelchair	33	Portable Suction Unit	7
Bath Bench	26	Power Adaptor	24
Bath Bench with Foldable Back Rest	27	Power Wheelchair	34
Bath Bench without Back Rest	27	Professional Stethoscope	22
Bathroom Weighing Scale	25	Professional Stethoscope (Cardio Model)	21
Bedding with Mattress	15	Professional Stethoscope (Sprague Rappaport)	21
Bedside Cabinet	16	Purifier Bedroom 402 Sandstone	20
Breathing Machine CPAP	11	Purifier Healthmate 400 Sandstone	20
Commode Chair without Wheels	27	Quad Cane	29
Commode with Wheels	27	Reclining High-Back Wheelchair	35
Cross Spoke Wheelchair	32	Reclining High-Back Wheelchair (with Commode)	35
Desk/Wall Type Sphygmomanometer	22	Replacement 5-Stage Filter for Bedroom	20
Defibrillator S5	11	Replacement Cover for Stethoscope Chestpiece	22
Digital Blood Pressure Monitor	23	Replacement Filter for Healthmate	20
Digital TENS Dual Channel	26	Rigid Walker	30
Double Step Foot Stool	17	Rubber Sheet	17
Dual Head Stethoscope	21	Rubber Tip for Folding Walker	30
Egg Cushion Mattress	15	ECG Machine	5
Emergency Carts	19	Service Cart	18
Examination Couch	15	Single Step Foot Stool	18
F9 Express Fetal & Maternal Monitor	6	Single-Use Tamper-Evident Lock	19
First Aid Kit 303	26	SonoTrax Basic A Fetal Doppler	6
Foldable Commode Chair	28	Spectacle Goggles	25
Foldable Walking Stick (Black)	29	Spine Board	18
Forearm Crutches	30	Stand-by Type Sphygmomanometer	23
Glass Breast Pump	26	Steel Wheelchair	33
Heavy Duty Suction Machine	8	Stretcher Patient Transfer Trolley	14
Heavy Duty Suction Unit	7	Suction Jar Lid Assembly	8
Heavy Duty Suction Unit	8	Syringe Pump	9
Ice Bag	24	TC-260 Oxygen Concentrator with Nebulization Function	12
iM50 Patient Monitor	6	Transit Aluminum Wheelchair	33
Infusion Pump	8	Ultrasonic Fetal Doppler	6
Infusion Pump	9	Underarm Crutches	30
IV Pole Stand	16	Walker	31
IV Stand for Wheelchairs (Two Prongs)	35	Walker with Wheels	31
Lightweight Mobility Scooter	35	Walking Stick (Offset Handle)	28
Mag Wheel Wheelchair	32	Walking Stick (Two Sections)	28
Manual Hospital Bed Side Railings Set	17	Walking Stick Flower Offset	29
Manual Up/Down Emergency Stretcher	13	Walking Stick Flower Offset, Two Sections	29
Medical Oxygen Regulator	11	Wheelchair Commode	33
Nebulizer Machine	10		
Nebulizer Machine	10		

MEDICAL EQUIPMENTS





ECG Machine 1-Channel

Edan Single Channel SE-1 is a compact single-channel ECG machine that record ECG waveforms in real-time, it displays clear traces with real-time heart rate and records real-time traces of user-selected leads.

FEATURES AND BENEFITS:

- ▶ Vivid waveform with foldable LCD display
- ▶ Built-in high resolution thermal printer
- ▶ Well-designed membrane keyboard
- ▶ Digital filters, resisting baseline drift, AC and EMG interference
- ▶ Automatic baseline adjustment
- ▶ Lead-off detection and alarm



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-EP-ECGMS1	Edan ECG Machine	Unit	1	1	MDRE




ECG Machine 3-Channel

The Edan ECG Machine is a high-quality and reliable three-channel ECG unit that features a compact and lightweight design, portable for any work situation.

FEATURES AND BENEFITS:

- ▶ Ultra-compact and lightweight design
- ▶ High resolution color touch screen
- ▶ 12-lead simultaneous acquisition and display
- ▶ Support linear/2D barcode scanner
- ▶ Signal quality identification
- ▶ Long time sampling
- ▶ Built-in rechargeable Li-ion battery



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-EP-ECGMS3	Edan ECG Machine	Unit	1	1	MDRE



ECG Machine 12-Channel

The Edan ECG Machine Express is designed to meet high-standard ECG requirements, designed with complete powerful data analysis functions, this unit offers reliable ECG solutions for both resting ECG and stress test.

FEATURES AND BENEFITS:

- ▶ Large high resolution colored touchscreen
- ▶ Comprehensive filters & anti-noise technology
- ▶ Real-time ECG waveforms freezing and review
- ▶ Advanced SEMIP algorithm for automatic measurements and interpretation
- ▶ Built-in high resolution thermal printer



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-EP-ECGM12	Edan ECG Machine	Unit	1	1	MDRE



Patient Monitor

Patient Monitor is a lightweight and portable patient monitor, a great choice for patient transport, with optional IBP and capnography support expands its applications to a wider usage in sub-acute divisions.

FEATURES AND BENEFITS:

- ▶ Optimized portable design for easy transportation
- ▶ Pacemaker detection, electrosurgical interference proof
- ▶ Defibrillation protection and defibrillation synchronization
- ▶ Unique iSEAP algorithm, optimized for arrhythmia patients
- ▶ Dual-mode anti-interference pulse oximeter
- ▶ Support barcode scanner
- ▶ Complete and flexible mounting solutions
- ▶ Nurse call function
- ▶ Standard Parameter: 3/5-lead ECG, HR, RESP, EDAN SpO2, NIBP, PR, 2-TEMP

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-EP-EPMM50	Patient Monitor	Unit	1	1	MDRE



Ultrasonic Fetal Doppler

TMS Ultrasonic Fetal Doppler is a high-performance fetal heartbeat detector, designed to meet the routine examination requirements of obstetricians.

FEATURES AND BENEFITS:

- ▶ Portable, lightweight and easy-to-use
- ▶ Full color OLED display FHR value, FHR intensity, improper FHR alarm signal, speaker volume, battery capacity indicator, working mode, frequency of the probe and FHR synchronization symbol
- ▶ Provides accurate FHR detection, it has FHR mode, Curve Mode, Replay Mode, Manual Count Mode, Average FHR Mode and Countdown FHR Mode
- ▶ With power saving that turns on automatically without signal within 30 seconds
- ▶ User-friendly machine interface for easy menu selection
- ▶ Standard 9V alkaline battery, available for power supply, which can support 300 minutes of continuous work

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TMS100	TMS Ultrasonic Fetal Doppler	Unit	1	20	—



Fetal Doppler

Fetal Doppler is a user-friendly digital heart monitoring device that includes auto and manual modes to select based on user's preference. The auto mode uses real-time FHR and average FHR measurement configurations for fetal heartbeat monitoring. With the manual mode, a view of the fetal heart rate updates in 5 or 10-second intervals for a calculated measurement.

FEATURES AND BENEFITS:

- ▶ Ergonomic design, compact and lightweight
- ▶ Supports 2MHz, 3MHz, 4MHz, 8MHz waterproof probes
- ▶ For automatic and manual FHR Counting
- ▶ Large backlight LCD display
- ▶ Ideal for long time and continuous use (8 hours)
- ▶ Accurate FHR detection with clear sound
- ▶ Automatic power-off
- ▶ Easy-to-use

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
P-MC-000006	Fetal Doppler	Unit	1	10	MDRE



Fetal & Maternal Monitor

Fetal & Maternal Monitor is an innovative fetal and maternal monitor that features integrated monitoring of twin's FHR, uterine activity, fetal movement, intrauterine pressure, direct ECG as well as maternal NIBP, SPO2, ECG and temperature. Designed for the obstetrical care continuum, spanning from the antepartum and intrapartum state to postpartum state.

FEATURES AND BENEFITS:

- ▶ 12.1" high resolution, colored TFT-LCD screen
- ▶ 12 crystal beam waterproof FHR transducer 12-crystal broad beam
- ▶ 24 hours back up memory for seamless monitoring
- ▶ Three monitoring interfaces for different clinical needs
- ▶ Signal overlap verification to separate twins' FHR
- ▶ Flat surface design with waterproof TOCO transducer

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-EP-EEFMF9	Fetal & Maternal Monitor	Unit	1	1	MDRE





Portable Defibrillator

TMS Portable Defibrillator is a compact, portable and lightweight defibrillator, designed with built-in module to ensure stable and reliable performance.

FEATURES AND BENEFITS:

- Multiple measuring functions which include: 3-lead, 7-lead, ECG/HR
- Complete built-in module design ensures stable and reliable performance
- Built-in recorder supports real-time recording, present screen printout and trigger print out by alarm
- Parameter display with big character
- 7-color high brightness TFT LCD monitor
- Portable design, stylish and convenient
- Rechargeable maintenance-free battery
- Resistant to high-frequency electrotome and is protected against defibrillation effects.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TD1000	TMS Portable Defibrillator	Unit	—	1	—



Heavy Duty Suction Unit

TMS Heavy Duty Suction Unit is a portable medical suctioning device that produces maximum vacuum of 560mmHg. The unit is equipped with 800ml collection bottle, designed for convenient and reliable usage.

FEATURES AND BENEFITS:

- New anti-bacterial filter design, helps prevent contamination
- High-quality vacuum regulator and anti-vibration vacuum gauge setting
- Oil-less permanent lubrication
- 800ml collection bottle with overflow protection

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TMS300	TMS Heavy Duty Suction Unit	Unit	1	1	MDRE



Portable Suction Machine

Thomas Portable Suction Machine is portable suctioning device that features high-performance, designed for applications requiring high vacuum capacity.

FEATURES AND BENEFITS:

- Oil-free diaphragm
- Capacity of plastic bottle cc/type: 1200cc/Autoclavable
- Hydrophobic filter to protect the pump
- Dimension: 37.5cm x 21cm x 25.5cm
- Max. vacuum: 560mmHg @ sea level
- Weight: 4.9 kg

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-TI-TP1637	Thomas Portable Suction Machine	Unit	1	1	MDRE



Portable Suction Unit

TMS Portable Suction Unit is a new generation oil-free suction device, specially designed to provide efficient vacuum operations.

FEATURES AND BENEFITS:

- Newly designed oil free lubrication pump
- Extremely portable for easy transport in the hospital, can also be used in home care settings
- With 1 piece 800cc collection bottle and 1 piece bacterial filter
- Sound level: 65dB
- Vacuum rate: 18LPM

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TMSSM1	TMS Portable Suction Unit	Unit	1	2	MDRE





Heavy Duty Suction Unit

TMS Heavy Duty Suction Unit is a mobile medical suctioning unit that produces maximum vacuum pressure of 630mmHG. The unit can work with two 3000ml collection bottles, four 45cm tubing and two 180cm tubing, and has been designed for convenient reliable usage and continuous operation.

FEATURES AND BENEFITS:

- ▶ New anti-bacterial filter design that helps prevent contamination
- ▶ High-quality vacuum regulator with anti-vibration vacuum gauge setting
- ▶ Collection bottle with overflow protection
- ▶ Quiet and effective operation due to the unique design of the WOBL Piston

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TMS520	TMS Heavy Duty Suction Unit	Unit	1	1	—



Heavy Duty Suction Machine

Thomas Heavy Duty Suction Machine Series is the perfect solution for different suction needs. This pump is a professional compact aspirator that is fully portable but can also be adapted as a table top unit.

FEATURES AND BENEFITS:

- ▶ Oil-free piston with a power of 1/7 HP
- ▶ Capacity of glass bottle cc/type: 2800cc/Autoclavable
- ▶ Hydrophobic filter to protect the pump
- ▶ Vacuum: 660mmHg
- ▶ Weight: 22.23kg
- ▶ Height: 90cm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-TI-HDS124	Thomas Heavy Duty Suction Machine	Unit	1	1	—



Suction Jar Lid Assembly

Thomas Suction Jar Lid Assembly are sturdy glass jar for suction pump with 2000cc capacity.

FEATURES AND BENEFITS:

- ▶ Autoclavable glass jars
- ▶ Capacity of 2000cc

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-TI-662984	Thomas Suction Jar Lid Assembly, 2000cc	Unit	1	1	—



Infusion Pump

TMS Infusion Pump is a reliable medical device that delivers fluids and medications into the patient's body in controlled amounts. Easy to operate by just turning the rotary dial to select desired settings. Compact and lightweight for easy handling and mounting.

FEATURES AND BENEFITS:

- ▶ Compatible with IV sets of any standard
- ▶ Real-time infused volume
- ▶ Automatic power switching
- ▶ Curvilinear peristaltic pumping mechanism
- ▶ With history log of 20,000 events
- ▶ With battery life of 5 hours at 25mL/hr
- ▶ With built-in thermostat (from 30-45°C)

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-TELF600	TMS Infusion Pump	Unit	1	1	MDRE





Infusion Pump

TMS Infusion Pump is an advanced infusion device, designed to guarantee accuracy of dose rate, engineered with; drug library, history record review and last configuration memory.

FEATURES AND BENEFITS:

- Universal: Compatible with any brand of IV sets
- Accuracy: High-precision performance
- Easy to use with intuitive operation system
- Safety: Alarm volume adjustable in 9 levels and bubble detection with 3 sensitivities
- Smart touch button, TFT oversize touchscreen
- With built-in thermostat (from 30-45 degree Celsius)

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-PM-TIP820	TMS Infusion Pump	Unit	1	1	MDRE



Syringe Pump

TMS Syringe Pump is an advance administration device that can accurately control syringe velocity, designed to administer fluid/medication to the patient in guaranteed accurate dose rate.

FEATURES AND BENEFITS:

- Universal: Compatible with any syringe brand, with different sets of standards
- Accuracy: People-oriented design, allows interactive operation, provides high-precision performance
- Easy-to-use: Intuitive operation system, drug library, history record review and last configuration memory
- Safety: Alarm volume adjustable in 9 levels, double CPU, bubble detection with 3 sensitivities
- Features: Smart touch button, high-definition touchscreen

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
4-PM-TSP900	TMS Syringe Pump	Unit	1	1	MDRE



Syringe Pump

KellyMed Syringe Pump is solidly built, this syringe pump can be relied on to infuse accurate doses of intravenous injectables over predefined time periods.

FEATURES AND BENEFITS:

- Applicable syringe: 12 brands for user-defined
- Occlusion sensitivity is adjustable
- Can be stacked up and integrated into multi-channel

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-KM-SP0602	KellyMed Syringe Pump	Unit	1	1	MDRE



Nebulizer with children inspired designs

Getwell Nebulizers with children inspired designs are shaped in fun, kid-friendly design for enjoyable and effective nebulization.

FEATURES AND BENEFITS:

- One button operation
- With overheat protection / auto shut-off
- Compact structure
- Portable and with carrying bag

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Bebe Penguin, Molly Cow & Jerry Fish	Unit	1	6	MDRE





Pocket Mesh Nebulizer

Getwell Pocket Mesh Nebulizer (Gen 1) is a small, portable and ultra-light nebulizer, designed to provide an effective and convenient nebulization treatment.

FEATURES AND BENEFITS:

- Quiet operation
- Lightweight and handy for effective nebulization anytime, anywhere.
- Provide efficient treatment, with low medicine residue.
- Compact design

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-IC-GWPMN0	Getwell Pocket Mesh Nebulizer	Unit	1	30	MDRE



Pocket Mesh Nebulizer

Getwell Pocket Mesh Nebulizer (Gen 2) offers the convenience of ultrasonic portables, allowing users to perform breathing treatments away from home. Produces the best quality aerosol mist that can be inhaled deeply into the respiratory tract for better control of respiratory ailments.

FEATURES AND BENEFITS:

- Lightweight and handy
- Allows fast delivery of medication
- Low noise operation
- Easy-to-clean

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-GWMN01	Getwell Pocket Mesh Nebulizer	Unit	1	30	MDRE



Nebulizer Machine

Getwell Nebulizer GWN-1250 is a compact nebulizer system with sleek and lightweight design that features a built-in nebulizer dock station that holds the nebulizer in place, preventing medication spillage.

FEATURES AND BENEFITS:

- Lightweight, compact design
- High-quality, oil-free, durable motor
- Silent, yet powerful piston operation
- With built-in handle for easy transport
- Easy-to-use with one button operation
- With built-in nebulizer dock that stabilize the nebulizer cup
- Ensures effective and fast medication delivery

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-GW1250	Getwell Nebulizer GWN-1250	Unit	1	8	—



Nebulizer Machine

Getwell Nebulizer GWN-1000 is a compact nebulizer system with sleek and lightweight design that features a built-in nebulizer dock station on one side, securing the nebulizer in place, preventing medication spillage.

FEATURES AND BENEFITS:

- One button operation
- With overheat protection/auto shut-off
- Compact structure
- Portable and with carrying bag

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-GWACN1	Getwell Nebulizer Machine GWN-1000	Unit	1	8	MDRE





Nebulizer Pump

TMS Nebulizer Pump is designed to deliver aerosol medication therapy effectively and efficiently with clinical accuracy. It is easy to use with its one-button operation and is ideal for all ages. The built-in handle allows easy transport from one place to another.

FEATURES AND BENEFITS:

- Effective medication delivery
- One-button operation
- Complete nebulizer kit included
- Easy-to-clean
- Lightweight

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-IC-TMSNEB	TMS Nebulizer Pump	Unit	1	8	MDRE



Breathing Machine CPAP

Yuwell Breathing Machine is a positive pressure breathing machine, used for treatment of adult sleep apnea syndrome.

FEATURES AND BENEFITS:

- 3.5" screen
- Automatic leakage compensation technology
- Pressure range: 4-20cmH2O
- PCM: noise removal technology
- Small and portable
- Automatic start and stop

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-YBM001	Yuwell Breathing Machine CPAP	Unit	1	1	—



Medical Oxygen Regulator

Inmed Medical Oxygen Regulator OXR-580 is a piston type flow meter regulator with easily accessible flow adjustment knob, designed to regulate the free flow of medical oxygen from an oxygen cylinder.

FEATURES AND BENEFITS:

- Piston type, high-pressure chamber with chromed forged brass body
- Polycarbonate inner and outer tubes
- Flow range: 0-15 LPM
- Working Pressure: 50PSI
- Back pressure compensated flow meter
- 3000 psi maximum inlet pressure
- Accessories: humidifier bottle, cannula and wrench
- Inlet Connection (connecting thread): CGA540
- Tested to CGA (Compressed Gas Association) standards

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-OXR580	Inmed Medical Oxygen Regulator	Unit	1	20	—



Portable Oxygen Concentrator (TC-170)

TMS Portable Oxygen Concentrator TC-170 is a lightweight, portable oxygen concentrator that generates oxygen and anion particles, kills bacteria and moisturizes air, creating pure and fresh air.

FEATURES AND BENEFITS:

- Stable oxygen concentration
- Delivers high oxygen concentration
- Ultra silent operation: 40dB
- Intelligent surveillance monitoring

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TMSPOC	TMS Portable Oxygen Concentrator TC-170	Unit	1	1	—





Oxygen Concentrator with Nebulization Function

The TMSOxygen Concentrator with Nebulization Function collects high purity oxygen from ambient air and is designed to provide assurance of oxygen supply.

FEATURES AND BENEFITS:

- Provides direct to patient high oxygen concentration.
- With fluctuation control feature, an innovative feature that allows oxygen output purity fluctuation rate controlled within 0.1%
- 6" LED screen display with E-touch operating panel
- Max oxygen flow: 5LPM, Oxygen Purity: 93%±3%
- Allows usage of up to 24 hours a day, 7 days a week
- Lightweight, movable with castors for convenience. Actual weight: 19kg
- 6 Safety Alarm System: low oxygen, oxygen pressure protection, pressure faulty, power off, compressor faulty, overheat protection

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
9-IC-TSM260	Oxygen Concentrator with Nebulization Function	Unit	1	—



Oxygen Concentrator (10L)

Angelbiss 10L Oxygen Concentrator is a high-volume oxygen concentrator that provides sufficient oxygen flow and efficient oxygen purity, designed for patients who require for more than 5LPM of oxygen flow rate. This unit features nebulization option and 6 alarm functions.

FEATURES AND BENEFITS:

- Touchscreen multi-function 6" LED light display
- With fluctuation control function that allows oxygen output purity fluctuation rate controlled within 0.1%
- High oxygen flow with oxygen concentration at 95%, all purity tested over 95%
- Single and dual flow options, two bottle system for a total of 10LPM
- Dual flow oxygen outlet accommodates two users at the same time
- Oxygen purity display by numbers with 96% accuracy
- 6 safety alarm function: Low Oxygen Alarm, Oxygen Pressure Protection Alarm, Pressure Faulty Alarm, Compressor Faulty Alarm, Overheat Protection and Power Off Alarm
- With nebulization function

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-TOC10L	Angelbiss Oxygen Concentrator, 10L	Unit	—	1	—



COMEN

Save lives with the reliable Comen Defibrillators

DEFIBRILLATOR S5

Engineered with comprehensive function that includes defibrillation, pacing, monitoring and AED mode. Features on-demand pacing and fixed pacing mode for patients in cardiac arrest or acute/severe arrhythmia. Pacing mode is rapid, easy-to-use, time-saving and improve recovery success rate.



Also available in:

S8 Model





2 Crank Manual Hospital Bed

2 Crank Manual Hospital Bed is a two control manual bed that perform back raise and leg raise functions. Height of the bed is fixed to accommodate most wheelchairs and commodes for easy patient transfer.

FEATURES AND BENEFITS:

- ▶ Manually adjustable backboard and foot board
- ▶ Allows quiet and smooth operation that provides maximum patient comfort
- ▶ Aluminum folding guard rails and pedal locking castors
- ▶ Two adjustment metal net board with detachable ABS head and end bed-boards
- ▶ Powder coating steel frame with double handle knobs



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-MB02A	Inmed 2 Crank Manual Hospital Bed	Unit	1	1	—



3 Crank Manual Hospital Bed with IV Stand

3 Crank Manual Hospital Bed is designed to provide safety and comfort for both patients and healthcare providers, it perform the following base functions: back raise, leg raise and height adjustment.

FEATURES AND BENEFITS:

- ▶ Manually adjustable backboard, calf board and foot board
- ▶ Quiet and smooth operation that provides maximum patient comfort
- ▶ Aluminum folding guard rails and pedal locking castors
- ▶ Two adjustment metal net board with detachable ABS head and end bed-boards
- ▶ Powder coating steel frame
- ▶ Double handle knobs
- ▶ With IV stand



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-MB03A	3 Crank Manual Hospital Bed with IV Stand	Unit	1	1	—

Stretcher Patient Transfer Trolley

Stretcher Patient Transfer Trolley is used for easy patient transfer and rescue, features an advanced centrally controlled brake system for stable and reliable operation.

FEATURES AND BENEFITS:

- ▶ High strength aluminum alloy material
- ▶ There are 2 height levels: fully elevated and ground level
- ▶ Adjustable backrest angle from 0° to 75°
- ▶ Product size: 185 (L) x 51 (W) x 24cm (H)
- ▶ Folded size: 92 (L) x 51 (W) x 11cm (H)

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-EQ-SPPT00	Stretcher (Patient Transfer Trolley)	Unit	1	1	—



3 Function Electric Hospital Bed with IV Stand

3 Function Electric Hospital Bed with IV Stand and Aluminum Alloy Guardrail provides both safety and comfort to patients and caregivers alike, whether it is in intensive care, standard operations and in patient transport.

FEATURES AND BENEFITS:

- ▶ Quiet and smooth operation that provides maximum patient comfort
- ▶ Electrically adjustable backboard, calf board and foot board
- ▶ Aluminum folding guard rails and pedal locking castors
- ▶ Collapsible aluminum siderail
- ▶ With double handle knobs
- ▶ With IV stand

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-OH-EHBIV3	3 Function Electric Hospital Bed with IV Stand	Unit	1	1	—



5 Function Electric Hospital Bed with IV Stand

5 Function Electric Hospital Bed with IV Stand provides the ultimate comfort level by allowing a multitude of adjustable positions which includes: back adjustment, leg adjustment, height adjustment, Trendelenburg and reverse Trendelenburg.

FEATURES AND BENEFITS:

- ▶ Overall length: 217cm
- ▶ Overall width: 98cm
- ▶ Lowest height: 35cm (floor to bed base)
- ▶ Highest height: 78cm (floor to bed base)
- ▶ Back raise angle: 0° to 70°
- ▶ Leg raise angle: 0° to 40°
- ▶ Tilt angle: 0° to 10°
- ▶ Mattress length: 200cm
- ▶ Mattress width: 90cm
- ▶ Safe working load: 160kg

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-OH-EHBIV5	5 Function Electric Hospital Bed with IV Stand	Unit	1	1	—



Manual Up/Down Emergency Stretcher (OH-TS014)

Manual Up/Down Emergency Stretcher OH-TS014 is used for patient transfer and rescue, with advanced centrally controlled brake system for stable and reliable operation.

FEATURES AND BENEFITS:

- ▶ Quality coated steel plate cart frame
- ▶ International advanced and centrally controlled brake system
- ▶ Adjustable back section with pneumatic spring
- ▶ Single crank for height adjustment

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-OH-MUDES0	Manual Up/Down Emergency Stretcher OH-TS014	Unit	1	1	—





Anti-Bedsore Mattress System Lattice Style

TMS Anti-Bedsore Mattress System Lattice style is made specially to prevent bedsore with its low pressure and large flow feature.

FEATURES AND BENEFITS:

- Provides excellent pressure reduction through continuous movement
- Clinically beneficial and suitable for both healthcare environment and home care settings
- Increases blood circulation
- Pneumatic type



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-PM-TMSADM	TMS Anti-Bedsore Mattress System, Lattice Style	Unit	1	1	—

Examination Couch Stainless Steel

Examination Couch Stainless Steel combines comfort, stability, strength and aesthetic design.

FEATURES AND BENEFITS:

- Sturdy and stable design
- Stainless steel frame
- Single hand operation
- Adjustable backrest
- Quick assemble design
- Product dimension: 190 x 60 x 65cm



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-EQ-EXMTRO	Examination Couch, Stainless Steel	Unit	1	1	—

Bedding with Mattress Uratex Foam

Bedding with Mattress is a pressure reduced foam with anti-microbial property.

FEATURES AND BENEFITS:

- Pressure reduced foam
- Low-friction and shearing force to protect the patient from skin damage
- Breathable water-resistant mattress covered with microbial inhibiting properties
- Easy cure covers, wipe clean and stain resistant




REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-MS-BMUF00	Bedding with Mattress, Uratex Foam	Unit	1	1	—

Egg Cushion Mattress

Egg Cushion Mattress is a latex-free cushion that provides comfort and pressure relief, prevents and helps heal skin ulcers, eliminates bed sores, reduces pressure pains, allows air to flow and prevents body from becoming too hot during sleep.

FEATURES AND BENEFITS:

- Pressure reduced foam
- Low-friction and shearing force to protect the patient from skin damage
- Breathable, water-resistant mattress cover with microbial inhibiting property
- Easy cure covers, wipe clean and stain resistant



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-MS-723016	Egg Cushion Mattress	Ea	1	1	—

INMED[®]

Bedside Cabinet

Inmed Bedside Cabinet is an economical and durable bedside cabinet, ideal for hospital rooms or home care use.

FEATURES AND BENEFITS:

- ABS plastic material
- Castor wheels
- Two layer shelves
- Two hanger
- Dimension: 48 x 48 x 76cm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-SBCL02	Inmed Bedside Cabinet	Unit	1	1	—

**INMED[®]**

IV Pole Stand

IV Pole Stand holds IV solutions securely, constructed from chrome plated steel, with adjustable height and 4 wheel base for stability.

FEATURES AND BENEFITS:

- Low profile base prevents tripping
- Allows transport of the pole in all surface
- Swivel ball bearing casters facilitates smooth roll
- Adjustable height: 47in. to 85in.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-IS4L4P	IV Pole Stand with 4 Legs & 4 Prongs	Ea	1	10	—
9-IC-IS4L2P	IV Pole Stand with 4 Legs & 2 Prongs	Ea	1	10	—

**INMED[®]**

Overbed Table

Inmed Overbed Table or hospital tray table allows people in bed to place food and their necessities within easy reach - right over them.

FEATURES AND BENEFITS:

- It has roll-on casters to move easily into position
- Allows patients to eat more comfortably and easily while in bed
- Finger tip control for height adjustment
- Spring loaded mechanism locks in position with slightest downward pressure

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-SCOT05	Inmed Overbed Table	Unit	1	1	—

**INMED[®]**

Overbed Table

Overbed table or hospital tray table allows people in bed to place food and their necessities within easy reach - right over them.

FEATURES AND BENEFITS:

- Light downward touch on table top locks table in position
- Tabletop is locked securely when height adjustment is released
- Easy to assemble with adjustable height
- Caster wheel for mobility

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-OH-OBTAB0	Overbed Table	Unit	1	2	—



Manual Hospital Bed Side Railings Set

Manual Hospital Bed Side Railings Set is designed to provide safety to patients by preventing risk of falling from the hospital bed, it also assists the patient in repositioning in the bed and in transitioning into or out of the bed.

FEATURES AND BENEFITS:

- Provides comfort and security
- Consists of 6 aluminum bars
- Approx. Dimension: 1440mm x 410mm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-0H-SROMHB	Manual Hospital Bed Side Railings Set	Pair	—	2	—



Rubber Sheet

Rubber Sheet is made from high-quality material; it features an exceptional high breaking and tearing strength.

FEATURES AND BENEFITS:

- Made from high-quality fabrics
- Durable and flexible
- Water-resistant

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-GE-RS2500	Rubber Sheet, 25 yards	Roll	1	—	—



Ordinary Castor

Ordinary Castor is a rubber caster which is versatile and can be utilized not only in professional environments, but also for home furnitures.

FEATURES AND BENEFITS:

- Quick and easy mounting
- Capacity of 350 pounds
- Heavy-duty solid rubber tires
- Durable 5 inches in height

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-0H-OCAST5	Ordinary Castor, 5 inches	Ea	1	—	—



Double Step Foot Stool

2 Layers

Double Step Foot Stool is handy and helpful in reaching high places. The first step is 10" high from the floor and the second step is 17 1/4" high from the floor.

FEATURES AND BENEFITS:

- Durable chrome finish
- Skid resistant rubber mat surface
- Extra-wide and deep tops
- Step platforms: 51x46x35cm
- Reinforced base and top

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-DSSTL0	Double Step Foot Stool, 2 Layers	Unit	1	1	—



Single Step Foot Stool

Single Step Foot Stool is a step-up stool, useful for assisting people, provides stability in reaching high places.

FEATURES AND BENEFITS:

- Extended offset legs, guards against tipping
- Deluxe model features chrome platform with rubber mat
- Rubber foot caps for stability
- Frame made from non-magnetic stainless steel pipes

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-EQ-SSSTL0	Single Step Foot Stool	Unit	1	1	—



INMED[®] Spine Board

Inmed Spine Board is a patient handling device, used primarily in pre-hospital trauma care. It is designed to provide rigid support during movement of a person with suspected spinal or limb injuries.

FEATURES AND BENEFITS:

- Constructed from high-density polyethylene
- Fully x-ray translucent
- Load bearing: 159kgs.
- Floatation capacity: 20kgs.
- Size: 185 x 45 x 6cm
- Weight: 8kgs.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-ISB001	Inmed Spine Board	Unit	1	1	—



Patient Lifter

Patient Lifter is an assistive device that allows transferring of patient with mobility problems from a bed or a chair into a similar resting place.

FEATURES AND BENEFITS:

- Power coated steel
- 150kg lifting capacity
- Long padded handles for various grip and for users of different height
- Directional wheel guide for straight steering, reduces stress to caregiver's back

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-EPL150	Patient Lifter	Unit	1	1	—



medbanc[®]

Service Cart

Medbanc Service Cart is a versatile and durable multi-purpose cart, designed to perform wide variety of tasks, with two shelves and one flat side for convenient loading.

FEATURES AND BENEFITS:

- Includes three shelves for storage
- Easy grip handles
- Easy to clean surfaces
- Measures: 15-1/8" x 19-1/8" x 29"
- Swivel casters for mobility

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-SCART0	Medbanc Service Cart	Unit	1	1	—





Emergency Carts with Tamper-Evident Lock

Medbanc Emergency Cart offers stable, smooth mobility and enhanced convenience, with changeable accessories, as well as tamper-evident lock for extra security.

FEATURES AND BENEFITS:

- Define a new standard of organization and function,
- With IV pole, defibrillator shelf, CPR board, oxygen tank holder and sharps disposal with bracket
- ABS top with slide out work surface
- Central lock control
- Intuitive design
- Available in blue and red

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-EC031B	Medbanc E-Cart 5 Drawers w/ Tamper-Evident Lock, 31 B	Unit	1	1	—
1-EQ-EC031R	Medbanc E-Cart 5 Drawers w/ Tamper-Evident Lock, 31 R	Unit	1	1	—
1-EQ-EC037B	Medbanc E-Cart 6 Drawers w/ Tamper-Evident Lock, 37 B	Unit	1	1	—
1-EQ-EC037R	Medbanc E-Cart 6 Drawers w/ Tamper-Evident Lock, 37 R	Unit	1	1	—
1-EQ-SACC00	Medbanc Standard Accessories for E-Cart	Set	1	1	—



Single-Use Tamper-Evident Lock

Medbanc Single-Use Tamper-Evident Lock is a tamper proof disposable lock, designed with one-time serial number, useful in emergency carts/medication kits for effective management of medication/emergency care items.

FEATURES AND BENEFITS:

- Material: PP (Polypropylene)
- Security tags: ideal for medication carts, crash carts, medication kit and other storage equipment that needs a safety lock.
- Each lock is printed with a unique serial number.
- Breakable without the use of any tools
- Sold in pack of 50's

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-MBSVPS	Medbanc Single-Use Tamper-Evident Lock	Pack	50	50	—



Purifier Bedroom 402 Sandstone

Austin Air Purifier Bedroom 402 Sandstone addresses the complete spectrum of air cleaning, removing sub-micron particles, noxious gases and chemicals. Austin Air's 360-degree intake system draws air into all sides of the the Bedroom Machine™, passing it through 5-stages of filtration for more cleaner air, delivered faster and efficiently.

FEATURES AND BENEFITS:

- ▶ A specialty air cleaner with 60sq. ft. of medical HEPA Filter, certified to remove 99.97% of airborne particles 0.3 microns or larger.
- ▶ Features 15lbs of granular carbon in addition to activated carbon cloth for better and more efficient absorption of odors and chemicals.
- ▶ Effectively cleans an area of 1500 sq. ft. (with normal ceiling heights).
- ▶ Effectively filters sub-micron particles, gases, odors and VOCs
- ▶ The included filter is designed to last up to five years in a residential setting.
- ▶ Comprehensive air purifier with five stages of filtration.
- ▶ Keeps your bedroom fresh and clean, developed to improve sleep.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-ARB402	Austin Air Purifier Bedroom 402 Sandstone	Unit	1	—	—



Purifier Healthmate 400 Sandstone

Austin Air Purifier Healthmate 400 Sandstone is a specialty air purifier with 60 sq. ft. of medical HEPA Filter (rated MERV-17), certified to remove 99.97% of airborne particles 0.3 microns and larger, ideal for improving indoor air quality especially in areas with poor air conditions.

FEATURES AND BENEFITS:

- ▶ Specialty air purifier with 60 sq. ft. of medical HEPA Filter (rated MERV-17), certified to remove 99.97% of airborne particles 0.3 microns and larger
- ▶ Contains 15 pounds of granular carbon and zeolite for the absorption of odors, gasses and chemicals found in smoke.
- ▶ Effectively cleans an area of 1500 sq. ft. (with normal ceiling heights).
- ▶ Cost economical. The included filter is designed to last up to five years in a residential setting.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-ARH400	Austin Air Purifier Healthmate 400 Sandstone	Unit	1	—	—



Replacement 5-Stage Filter for Bedroom 402

Replacement filter for the Austin Air Purifier Bedroom 402, features 5-stage filtration system for overall particulate protection.

FEATURES AND BENEFITS:

- ▶ STAGE 1 – Large Particle Pre-filter. Removes particles easily seen by the naked eye
- ▶ STAGE 2 – Medium Particle Pre-filter. Removes small to medium size particles
- ▶ STAGE 3 – Over 780 cubic inches of Activated Carbon and Zeolite Removes chemicals, gases and odors
- ▶ STAGE 4 – 60 sq. ft. of True Medical Grade HEPA Filter. Removes 99.97% of all particles larger than 0.3 microns.
- ▶ STAGE 5 – Military Grade HEGA Carbon Cloth. HEGA (High Efficiency Gas Arrestance), providing an extremely large surface area for adsorption of chemicals, gases, viruses and bacteria

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-ARFW05	Austin Air Replacement 5-Stage Filter for Bedroom 402	Ea	1	—	—



Replacement Filter for Healthmate 400

Replacement filter for the Austin Air Healthmate, constructed using a true medical High Efficiency Particulate Air (HEPA) filter (rated MERV-17) surrounded by a carbon/zeolite mixture, capable of removing the tiniest of particles and absorbing odors and gases. It uses HEPA to remove up to 99.97 percent of all particles down to 0.3 microns in size.

FEATURES AND BENEFITS:

- ▶ STAGE 1 – Large Particle Pre-filter. Removes particles easily seen by the naked eye
- ▶ STAGE 2 – Medium Particle Pre-filter. Removes small to medium size particles
- ▶ STAGE 3 – Over 780 cubic inches of Activated Carbon and Zeolite. Removes chemicals, gases and odors.
- ▶ STAGE 4 – 60 sq. ft. of True Medical Grade HEPA Filter. Removes 99.97% of all particles larger than 0.3 microns.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-ARFW00	Austin Air Replacement Filter for Healthmate 400	Ea	1	—	—



OTHER HEALTH CARE PRODUCTS

INMED[®]

Dual Head Stethoscope

Inmed Dual Head Stethoscope features a double sided chestpiece, engineered to provide maximum sound pick-up and transmission.

FEATURES AND BENEFITS:

- With eartips that fit comfortably
- Long lasting membrane on diaphragm
- Overall length is 28.5 inches
- Chrome-plated brass binaural tubing

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-IDHSTC	Inmed Dual Head Stethoscope	Ea	1	100	—



Brielle[®]

Professional Stethoscope (Sprague Rappaport)

Brielle Select Professional Stethoscope Sprague Rappaport offers improved acoustic performance. Single piece with two tubes that are connected by a continuous panel, this prevents the tubes from rubbing each other in use and making distracting noises, and it provides the stethoscope much more compact audibility of high-frequency sounds.

FEATURES AND BENEFITS:

- High-acoustic sensitivity
- Combination of stainless steel chestpiece with non-chill rim
- Stainless steel internal spring binaural
- High-quality extra thick Y PVC tubing
- With free spare ear tips and diaphragm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-LS2205	Brielle Select Professional Stethoscope, Sprague Rappaport	Ea	1	50	—



Brielle[®]

Professional Stethoscope (Cardio Model)

Brielle Select III Professional Stethoscope Cardio Model is used to hear sounds made by the heart, lungs, and intestines, offering improved acoustic performance with better audibility of high-frequency sounds.

FEATURES AND BENEFITS:

- High-acoustic sensitivity
- Combination of stainless steel chest piece with non chill rim
- Stainless steel internal spring binaural
- High-quality extra thick Y PVC tubing
- With free spare ear tips and diaphragm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-LS2401	Brielle Select III Professional Stethoscope, Cardio Model	Ea	1	25	—





Professional Stethoscope

Brielle Professional Stethoscope Select Model provides reliable acoustic performance with better audibility of high-frequency sounds.

FEATURES AND BENEFITS:

- High-acoustic sensitivity
- Combination of stainless steel chestpiece with non chill rim
- Stainless steel internal spring binaural
- High-quality extra thick Y PVC tubing
- With free spare ear tips and diaphragm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-SS2301	Brielle Select Professional Stethoscope, Adult	Ea	1	25	—
9-IC-LS2302	Brielle Select Professional Stethoscope, Pedia	Ea	1	25	—



Replacement Cover for Stethoscope Chestpiece

3M Littmann Replacement Cover for Stethoscope Chestpiece is a replacement diaphragm cover for Littmann stethoscope. Fits securely over stethoscope head. Replace regularly to prevent infection and to ensure high quality sound transmission.

FEATURES AND BENEFITS:

- Clear material that provides an airtight seal to the chestpiece.
- Diameter: 43mm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-MS-723001	3M Littmann Replacement Cover for Stethoscope Chestpiece	Ea	1	—	—



Aneroid Sphygmomanometer

Inmed Aneroid Sphygmomanometer type with dial reading up to 300mmHg is designed to provide accurate blood pressure monitoring, anytime, anywhere.

FEATURES AND BENEFITS:

- 300mmHg aneroid meter
- Velcro cuff
- Durable air chamber
- Pump and control

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-IASMMM	Inmed Aneroid Sphygmomanometer	Ea	1	—	MDRE



Desk/Wall Type Sphygmomanometer

Inmed Desk/Wall Type Sphygmomanometer allows convenient blood pressure monitoring with accurate result.

FEATURES AND BENEFITS:

- Desktop model
- May also be installed on the wall via screw (2-in-1 function)
- 14cm diameter dial
- 300mmHg aneroid meter
- Vertical angulation for easy reading
- Cotton cuff with Velcro
- Durable air chamber, pump and control

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-GC0502	Inmed Desk/Wall Type Sphygmomanometer	Ea	1	10	—



INMED[®]

Stand-by Type Sphygmomanometer

Inmed Stand-by Type Sphygmomanometer is a monitoring device with basket and castor, designed for convenient BP monitoring, allowing easy transport from one place to another.

FEATURES AND BENEFITS:

- Adjustable polycarbonate stand
- 14cm diameter dial
- 300mmHg aneroid meter
- Vertical angulation for easy reading
- Cotton cuff with velcro
- Durable air chamber, pump and control
- 5 wheels with brake
- Metal basket

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-GC0701	Inmed Stand-by Type Sphygmomanometer	Set	1	1	—

**Brielle[®]**

Automatic Arm-Type Blood Pressure Monitor (BM-26U)

Brielle Automatic Arm-Type Blood Pressure Monitor BM-26U is an arm type BP Apparatus, designed to provide convenient usage and accurate measurements.

FEATURES AND BENEFITS:

- Easy-to-use with large LCD display, measures blood pressure and pulse heart rate
- Large LCD screen display function, digital LCD display clearly shows systolic blood pressure, diastolic blood pressure, pulse rate, date, time and user measurements
- Ideal for two users
- Auto shut off feature
- Blood pressure meter using intelligent pressure automatic measurement technology
- High accuracy measurement with ± 3 mmHg for Pressure and $\pm 5\%$ for Pulse Rate
- Dual user mode, intelligent pressurization, two people can use, switch at any time
- Compact and lightweight
- With voice prompt function



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-BBM26U	Brielle Automatic Arm-Type Blood Pressure Monitor, BM-26U	Unit	1	20	MDRE

**Brielle[®]**

Automatic Arm-Type Blood Pressure Monitor (BM-63U)

Brielle Automatic Arm-Type Blood Pressure Monitor BM-63U is an arm type BP Apparatus, designed to provide convenient usage and accurate measurements.

FEATURES AND BENEFITS:

- Lightweight, handy and compact design
- LCD size is just 3.45" X 2.05"
- One button operation
- Silent inflation for convenient and discreet measurement
- Measurement range: 0-280mmHg (0-37.30kPa), Pulse: 40-199 beats/min
- Cuff Pressure Range: 0-299mmHg (0-39.86kPa)
- High accuracy measurement with ± 3 mmHg for Pressure and $\pm 5\%$ for Pulse Rate
- Dimension: 5.51"(W) X 2.89"(L) X 1.23"(H)



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-BBM63U	Automatic Arm-Type Blood Pressure Monitor, BM-63U	Unit	1	24	MDRE



Automatic Arm-Type Blood Pressure Monitor (BM-270)

Brielle Automatic Arm-Type Blood Pressure Monitor with visual color indicator for BP range, with time and date display for 2 person usage.

FEATURES AND BENEFITS:

- Large-scale LCD
- 22-45cm testing cuff
- Integrated air plug design
- Cuff wrapping detection
- IHB indication
- Body movement detection
- 74 groups of memory
- Average value display of three times measurement
- Adaptor input: AC 100-240V-50/60Hz 0.15A
- Adaptor output: DC 6V=60mA

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-BYE670	BM-270	Unit	1	24	MDRE



Automatic Arm-Type Blood Pressure Monitor (BM-130)

Brielle Automatic Arm-Type Blood Pressure Monitor with visual color indicator for BP range, time and date display for 2 person usage.

FEATURES AND BENEFITS:

- 75x54mm LCD display
- Detects cardiac anomalies
- Convenient for 2 users, stores up to 90 memory each
- Core bio-sensor technology to ensure accurate measurement
- Colored backlight to easily identify the range/level of the blood pressure result
- Averages the last three measurements

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-BAB806	BM-130	Unit	1	24	MDRE



Power Adaptor for Brielle BM-130

Power Adaptor compatible with Brielle digital blood pressure monitor.

FEATURES AND BENEFITS:

- Adaptor for Brielle digital blood pressure monitor
- AC 100-240V, 50/60Hz, 0.3A input
- DC 6V, 1A output

REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-BABPA1	Power Adaptor for Brielle BM-130	Unit	1	24	MDRE



Ice Bag

Inmed Ice Bag provides the right size to cover a decent size area for relieving swelling or muscular pain.

FEATURES AND BENEFITS:

- Pleats for expanded capacity
- Molds to fit body contour
- Extra large openings for easy refilling
- Available sizes: 6 and 9

DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
Inmed Ice Bag, 6 & 9	Ea	1	100	MDRE



INMED[®]

Spectacle Goggles

Inmed Spectacle Goggles provides eye protection against flying debris or against visible and near-visible light or radiation.

FEATURES AND BENEFITS:

- Built-in side shields and brow guard
- Excellent clarity
- Lightweight



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-ISBBT1	Inmed Spectacle Goggles	Ea	10	250	—

Bathroom Weighing Scale

Bathroom Weighing Scale is calibrated to measure body weight or calculated body mass.

FEATURES AND BENEFITS:

- High-quality ABS plastic
- Max weight: 130kg/280 lbs



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
0-GE-BR9010	Bathroom Weighing Scale	Unit	1	—	—

INMED[®]

Physician Weighing Scale

Height and Weighing Scale is the elegant all around solution for height and weight identification, suitable for use in the medical sector for doctor's clinics and hospitals.

FEATURES AND BENEFITS:

- High-quality ABS plastic
- Max weight: 130kg/280 lbs



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
6-GE-UPWS001	Height and Weighing Scale	Unit	1	—	—

INMED[®]

Baby Weighing Scale

Inmed Baby Weighing Scale provides security and support during infant weighing and scaling.

FEATURES AND BENEFITS:

- Maximum weight capacity: 20kg. / 44lbs.
- Plate dimension: 54.5 x 26cm
- Division: 50g (0~10kg) / 2oz(0~22lb)
- Overall dimension: 54.5 x 37 x 15cm
- 100g (10~20kg) / 4oz (22~44lbs)
- Net weight: 3.1kg.



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
9-IC-GC0002	Inmed Baby Weighing Scale	Unit	1	—	—



Digital TENS Dual Channel

Brielle Digital TENS Dual Channel is an electrical stimulator, that aims to provide a degree of symptomatic pain relief by exciting the sensory nerves.

FEATURES AND BENEFITS:

- ▶ Dual channel operation with six different stimulation modes and intensity
- ▶ Compact for better portability
- ▶ LCD display and digital operations
- ▶ Large buttons for easy operation
- ▶ Max. output 200Hz for more severe pain
- ▶ Intensity is adjustable to level 9



REF CODE	DESCRIPTION / SIZE	UOM	QTY/BOX	QTY/CTN	MDR
1-BR-BD7900	Brielle Digital TENS, Dual Channel	Unit	1	100	MDRE

Glass Breast Pump

Glass Breast Pump is a device used to extract milk from the breast of a lactating woman.

FEATURES AND BENEFITS:

- ▶ Maximize milk production and personal comfort level
- ▶ Has glass reservoir with rubber squeeze bulb



REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-MS-BPG000	Glass Breast Pump	Ea	1	—	—



First Aid Kit 303

TMS First Aid Kit with 303 essential first aid items. Contents are organized and labeled with clear directions for finding the right item for the injury. Contains concise directions for use.

FEATURES AND BENEFITS:

- ▶ Carrying case: Hard Case and Soft Case bag available
- ▶ 303 pieces, organized via the Easy Access Pocket System that allows faster and ease-of-use
- ▶ Components: Small Cuts and Burns, Medium Cuts and Scratches, Severe Bleeding and Burns, CPR, Protection and Instruments
- ▶ Clearly labeled pockets

Hard Case

Soft Case



DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Hard Case & Soft Case	Ea	1	10	—

Bath Bench with Back Rest

Bath Bench with Back Rest is made of blow molded high-density polyethylene plastic. Features backrest, drainage holes and built-in patient assist handles.

FEATURES AND BENEFITS:

- ▶ Seat and backrest made of blow molded high-density polyethylene plastic
- ▶ Drainage hole and built-in patient assist handles
- ▶ Adjustable height legs with outward angles and wide base to prevent accidental overturning
- ▶ Heavy duty anti-rust legs made from anodized aluminum tubing
- ▶ Extra large rubber tips for better safety and stability



REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-FS798L	Bath Bench with Back Rest	Ea	1	2	—

Bath Bench without Back Rest

Bath Bench without Back Rest is suitable for individuals who have no trouble with balance or sitting endurance.

FEATURES AND BENEFITS:

- Seat is made of blow molded high-density polyethylene plastic
- Drainage hole and built-in patient assist handles
- Adjustable height legs with outward angles and wide base to prevent accidental overturning
- Heavy duty anti-rust legs made from anodized aluminum tubing
- Extra large rubber tips for better safety and stability

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-FS797L	Bath Bench without Back Rest	Ea	1	2	—



Bath Bench with Foldable Back Rest

Bath Bench with Foldable Back Rest provides stability for the user for personal hygiene tasks like bathing. Bath bench enables individuals to be more independent without compromising safety.

FEATURES AND BENEFITS:

- Excellent quality blow molded high-density plastic seat for safety and comfort
- Anodized aluminum frame, angle front and rear legs that add stability
- Built-in patient assist handles with drainage holes
- Bench can be folded when not in use, easy to carry
- Big reinforced rubber tips, designed for stability

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-FS7900	Bath Bench with Foldable Back Rest	Ea	1	2	—



INMED[®] Commode with Wheels

Wheeled commode is designed for people with limited mobility. Used to transfer patient to and from any locations, with swivel casters allowing the commode to move effortlessly. Includes a white plastic seat with removable bucket and lid. Hygienic and easy to clean.

FEATURES AND BENEFITS:

- Easy detachable armrest and footrest for convenient patient transfer
- Powder coating wheelchair
- Standard plastic toilet seat, bucket and cover
- Detachable upholstered seat
- 2 swivel casters (with lockable brakes)

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-CW4LWF	Commode with Wheels	Ea	1	1	—



INMED[®] Commode Chair without Wheels

Inmed Commode Chair is designed for people with limited mobility. Equipped with an extra-large commode bucket that can be removed and emptied. Each leg has a non-slip tip to provide extra stability and safety. Comfortable and hygienic to use.

FEATURES AND BENEFITS:

- Powder coated chair
- Padded backrest and armrest
- Detachable upholstered seat
- Standard plastic toilet seat, bucket and cover

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-CWW893	Inmed Commode Chair Without Wheels	Ea	—	1	—



INMED[®]

Foldable Commode Chair

Inmed Foldable Commode Chair is designed for people with limited mobility. Equipped with commode bucket that can be removed and emptied. Each leg has a non-slip tip to provide extra stability and safety, with foam cushion for comfortable usage.

FEATURES AND BENEFITS:

- Powder coated chair
- Padded armrest
- Detachable upholstered seat
- Standard plastic toilet seat, bucket and cover

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-CR6621	Inmed Foldable Commode Chair w/ Foam & Bucket	Ea	—	1	MDRE

**INMED[®]**

Pedal Exerciser

Inmed Pedal Exerciser is used for upper and lower body exercise, improving circulation, muscle strength, joint range of motion and coordination.

FEATURES AND BENEFITS:

- Plastic pedals
- Chrome plated steel frame construction
- Adjustable resistant knob
- Comes fully-assembled

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-PE960F	Inmed Pedal Exerciser	Ea	1	—	—

**INMED[®]**

Walking Stick (Offset Handle)

Inmed Walking Stick with Offset Handle is used primarily to aid in walking, providing postural stability, balancing and support to the user.

FEATURES AND BENEFITS:

- Offset handle with foam and strap
- Adjustable height
- Chrome plated

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-WSOHD	Inmed Walking Stick Offset Handle	Ea	1	20	—



Walking Stick (Two Sections)

Walking Stick Two Sections is used primarily to aid in walking, providing postural stability, balancing and support to the user.

FEATURES AND BENEFITS:

- Offset handle with foam and strap
- Chrome plated
- Adjustable height

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-WSB250	Walking Stick, 2 Sections	Ea	1	—	—





Foldable Walking Stick (Black)

Foldable Walking Stick is used to support when walking.

FEATURES AND BENEFITS:

- ▶ Made of aluminum tubing
- ▶ High-quality hard-anodized treatment for avoiding scrape and surface damage
- ▶ Reinforced rubber tips
- ▶ Adjustable height
- ▶ Left or right made quickly by depressing locking button
- ▶ Folds away neatly for easy storage when not in use
- ▶ Max user weight 100kg.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-FS927L	Foldable Walking Stick, Black	Ea	1	20	—



Walking Stick Flower Offset

Walking Stick Flower Offset is a walking aid with flower design and an offset handle, designed to distribute the user's weight fully and evenly along the shaft of the cane all the way down to the tipped base.

FEATURES AND BENEFITS:

- ▶ Offset handle with foam and strap
- ▶ Chrome plated
- ▶ Adjustable height
- ▶ Flower Design

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-WSF000	Walking Stick, Flower Offset	Ea	1	—	—



Walking Stick Flower Offset, Two Sections

Walking Stick Flower Offset, Two Section is used to facilitate balancing while walking.

FEATURES AND BENEFITS:

- ▶ Offset handle with foam and strap
- ▶ Adjustable 80 - 90cm height
- ▶ Chrome plated
- ▶ Flower Design

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-WSF250	Walking Stick, Flower Offset Two Sections	Ea	1	—	—



Quad Cane

Quad Cane gives very firm and secure support to the user, it has four pronged base for increased stability and balance.

FEATURES AND BENEFITS:

- ▶ Made of aluminum material
- ▶ Contoured handle
- ▶ Fitted to left or right handed person
- ▶ Adjustable height 65 - 87.5cm
- ▶ Base length is 8"
- ▶ Base width is 6"
- ▶ Can support up to 850kg.

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-LBQCCF	Quad Cane	Ea	1	10	—



INMED

Forearm Crutches

Forearm Crutches is used to counter a mobility impairment or an injury that limits walking ability.

FEATURES AND BENEFITS:

- The ultimate lightweight forearm crutches that combines high-quality and good looks
- Contoured arm cuffs increase comfort and stability for the user, along with the general forearm crutch design

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Forearm Crutches, Adult & Child	Ea	1	20	—

**INMED**

Underarm Crutches

Underarm Crutches is a support fitting under the armpit for use as an aid in walking.

FEATURES AND BENEFITS:

- Aluminum material
- Comes in child, young adult and adult designs
- Solid and adjustable

DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
Underarm Crutches, Young & Adult	Ea	1	20	—



Rubber Tip for Folding Walker

Rubber Tip for Folding Walker offers extra reinforcement and stability for the walker.

FEATURES AND BENEFITS:

- Slip-resistant
- Shock Absorbing

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
0-MS-C06002	Rubber Tip for Folding Walker	Ea	1	—	—

**INMED**

Rigid Walker

Rigid Walker is a tool for disabled or elderly people who need additional support to maintain balance or stability while walking, features an ergonomic handle for comfortable grip.

FEATURES AND BENEFITS:

- Standard walking frame
- Cushioned handgrips
- 4 rubber tips for added comfort and stability

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-WREMFD	Rigid Walker, Medium	Ea	1	10	—
9-IC-WRELF	Rigid Walker, Large	Ea	1	10	—



INMED[®]

Walker

Walking Aid provides support to maintain balance or stability during walking.

FEATURES AND BENEFITS:

- Made of 7/8" lightweight anodized aluminum
- Has cushioned handgrips
- Dual purpose can be fixed or used as a reciprocal walker
- The large single button release can be operated even by those with limited dexterity
- Double cross brace for additional stability

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-S913L0	Walker	Ea	1	2	—

**INMED[®]**

Walker with Wheels

Inmed Walker with Wheels is specially designed for elderly person to provide walking support for stability and mobility. It has rubber grips for firm support and can be folded easily.

FEATURES AND BENEFITS:

- Alloy frame with ultra lightweight aluminum material
- Height adjustable: 30" to 37"
- Height adjustable in the leg area
- For folding action, push button is provided
- Slip-resistant

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-WWW912	Inmed Walker with Wheels	Ea	1	4	—

**INMED[®]**

2-Wheels Rollator

Rollator 2 Wheels is a rehabilitation rollator which gives recovering patients or aged people support and security to walk the distance without relying on other person.

FEATURES AND BENEFITS:

- Lightweight & durable
- With 2 rubber tips
- Foam padded backrest
- Adjustable height
- Handle with foam

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-2WR914	Rollator, 2 Wheels	Ea	1	2	—

**INMED[®]**

4-Wheels Rollator

Rollator 4 Wheels is a stylish and modern rolling walker that features a durable steel frame and comfortable seat.

FEATURES AND BENEFITS:

- Lightweight & durable
- Folds easily with basket attached
- Foam padded backrest
- With rubber tips
- Hand breaks on the handles for safety

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-4WR965	Rollator, 4 Wheels	Ea	1	1	—



INMED[®]

Pediatric Wheelchair

Inmed Pediatric Wheelchair is a portable pediatric wheelchair with perfect unsurpassed quality and style.

FEATURES AND BENEFITS:

- ▶ Aluminum lightweight
- ▶ Pneumatic rear wheels
- ▶ Fixed armrest and footrest
- ▶ Drop back handle
- ▶ Solid castors
- ▶ United brake

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-PWFS35	Inmed Pediatric Wheelchair	Ea	1	1	—



Pediatric Wheelchair

Pediatric Wheelchair is a portable pediatric wheelchair with perfect unsurpassed quality and style.

FEATURES AND BENEFITS:

- ▶ Children type aluminum chair frame
- ▶ Solid magwheels
- ▶ Fixed armrest and footrest
- ▶ Drop back handle
- ▶ Solid castors
- ▶ United brake

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-FS980L	Pediatric Wheelchair	Ea	1	1	—

**INMED[®]**

Cross Spoke Wheelchair

Inmed Cross Spoke Wheelchair is a device allowing manual propulsion by the seated occupant turning the rear wheel by hand. Cross spoke pattern makes the wheel strong and meets moderate to high durability.

FEATURES AND BENEFITS:

- ▶ Chrome steel frame
- ▶ Solid castors
- ▶ Fixed armrest and footrest
- ▶ Metal footrest with leg rest strap
- ▶ Foamed armrest
- ▶ Solid rear wheel with spokes

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-W80946	Inmed Cross Spoke Wheelchair	Ea	1	—	—

**INMED[®]**

Mag Wheel Wheelchair

Inmed Mag Wheel Wheelchair is a standard wheelchair with foldable chrome steel frame, designed for basic functionality, with solid 8 spoke magwheel composite construction for strength and durability.

FEATURES AND BENEFITS:

- ▶ Foldable chrome plated wheelchair
- ▶ Metal footrests with legrest straps
- ▶ Fixed footrest, armrest and backrest
- ▶ 18 inches seat width
- ▶ Solid mag wheel and castors
- ▶ 8-Spoke mag wheel

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-WC18FK	Inmed Mag Wheel Wheelchair	Ea	1	—	—



Steel Wheelchair

Steel Wheelchair is a stylishly designed wheelchair that ensures user's mobility and comfort.

FEATURES AND BENEFITS:

- Green powder coating steel frame
- Flip-up padded desk-length armrest
- Solid mag wheel
- Detachable footrest
- Solid castors

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
1-EQ-FS909B	Steel Wheelchair, 46cm/18" width	Ea	1	—	—



INMED[®]

Bariatric Wheelchair

Inmed Bariatric Wheelchair is a specially designed model with wide seat in 22" for bariatric user. It comes with durable powder coated steel frame that looks fashionable. The reliable wheelchair with dual cross brace offers a loading weight capacity that can carry up to 250lbs.

FEATURES AND BENEFITS:

- Flip up armrest, detachable legrest/footrest
- Powder coating steel frame
- PU Castor
- Seat width 22"
- Mag wheels solid tires

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-FWB951	Inmed Bariatric Wheelchair, 22"	Ea	1	—	—



INMED[®]

Wheelchair Commode

Wheelchair Commode FS609 is a standard wheelchair with commode that gives the user the independence on their mobility and other daily routines.

FEATURES AND BENEFITS:

- Foldable stainless steel frame
- Solid front castor wheel
- Fixed footrest and armrest
- Solid rear wheel
- 300lbs. weight capacity
- Plastic commode seat

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-WC609F	Wheelchair Commode FS609	Ea	1	1	—



INMED[®]

Transit Aluminum Wheelchair

Inmed Transit Aluminum Wheelchair is a light aluminum travel and transport wheelchair which features durable castors for propulsion. Economical and configured with fold feature for the easiest transportation.

FEATURES AND BENEFITS:

- Lightweight aluminum folding wheelchair
- United brake and seatbelt
- Flip-up armrest
- 6" Solid mag rear wheels
- Drop-back handle
- 5" Front castor wheels

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	QTY/CTN	MDR
9-IC-TAW001	Inmed Transit Aluminum Wheelchair	Ea	1	1	—



Aluminum Wheelchair

Aluminum Wheelchair is specially designed to be compact and travel-friendly. This model is best suited for attendant assisted especially for small adults and children.

FEATURES AND BENEFITS:

- Super lightweight / powder coating
- Fixed leg rest and armrest
- Solid castors
- Drop back handle with united break
- Pneumatic wheels
- Tubeless tires

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
1-EQ-FS871L	Aluminum Wheelchair 41cm/16" width, Child	Ea	1	—
1-EQ-FS871B	Aluminum Wheelchair 41cm/16" width, Adult	Ea	1	—



Aluminum Wheelchair (16" width)

Aluminum Wheelchair is specially designed to be compact and travel-friendly. This model is best suited for attendant assisted, especially for small adults and children.

FEATURES AND BENEFITS:

- Super lightweight / powder coating
- Fixed leg rest and armrest
- Solid castors
- Drop back handle with united break
- Pneumatic wheels
- Tubeless tires

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
1-EQ-S87446	Aluminum Wheelchair 41cm/16" width	Ea	1	—



Aluminum Wheelchair (18" width)

Aluminum Wheelchair is specially designed to be compact and travel-friendly. Suitable for hospital use and travels.

FEATURES AND BENEFITS:

- Quick release PU magwheel with safety and anti-tilt properties
- Detachable supporting back wheel
- Flip-up high adjustment desk armrest
- Detachable strong footrest
- Folding aluminum alloy
- Aluminum chair frame
- Solid magwheels
- With extra foam for comfort
- Solid castors

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
1-EQ-FS950L	Aluminum Wheelchair, 46cm/18" width	Ea	1	—



INMED[®]

Power Wheelchair

Inmed Power Wheelchair is a lightweight standard power wheelchair with easy to operate controller. Equipped with high powered motor batteries and intelligent electromagnetic brake which guarantees user safety.

FEATURES AND BENEFITS:

- Detachable armrest, footrest and rubber foam castor
- Powder coating steel frame
- Drop back handle
- Foldable chair frame
- Solid drive wheel with hand rim

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CTN	MDR
9-IC-PW2146	Inmed Power Wheelchair	Ea	1	—





Reclining High-Back Wheelchair (with Commode)

Inmed Reclining High-Back Wheelchair with Commode is specially designed to allow a person to comfortably recline at incremental angles: safely and comfortably. Equipped with an extra-large commode bucket that can be removed and emptied.

FEATURES AND BENEFITS:

- Foldable chromed steel frame
- PVC caster wheel and 24" solid wheel
- U-shaped commode seat
- Removable plastic commode pail with lid
- Detachable armrest, detachable elevating footrest
- Push to lock wheel brakes
- Reclining high backrest

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CARTON	MDR
9-IC-WRHBWC	Inmed Reclining High-Back Wheelchair w/ Commode	Ea	1	—



Reclining High-Back Wheelchair

Reclining High-Back Wheelchair is specially designed to allow a person to comfortably recline at incremental angles, safely and comfortably.

FEATURES AND BENEFITS:

- Foldable Chromed steel frame
- With mag wheels and 24" solid wheel
- Detachable armrest, detachable elevating footrest
- Push to lock wheel brakes
- Reclining high backrest

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CARTON	MDR
9-IC-WRHBFD	Reclining High-Back Wheelchair	Ea	1	—



Lightweight Mobility Scooter

Luggie Lightweight Mobility Scooter is a lightweight, portable and foldable mobility aid scooter, designed to increase independence and freedom for people who have poor walking and long-term health conditions.

FEATURES AND BENEFITS:

- Innovative and lightweight design
- Sturdy alloy aluminum chassis
- Folded scooter conveniently rolls along like a travel suitcase
- Adjustable tiller height and angle
- Foldable in 3 easy steps
- With 8.5Ah lithium battery

REF CODE	DESCRIPTION / SIZE	UOM	QTY/CARTON	MDR
1-EQ-LUG3A0	Luggie Lightweight Mobility Scooter, Black	Unit	1	—
1-EQ-LUG30R	Luggie Lightweight Mobility Scooter, Burgundy Red	Unit	1	—
1-EQ-LUG30Y	Luggie Lightweight Mobility Scooter, Yellow	Unit	1	—
1-EQ-LUG30C	Luggie Lightweight Mobility Scooter, Champagne	Unit	1	—



IV Stand for Wheelchairs (Two Prongs)

IV Stand for Wheelchairs with Two Prongs attaches to most wheelchairs, easily hangs IV bags/blood bag.

FEATURES AND BENEFITS:

- With height adjustment
- Chrome plated finish
- Attaches easily to any wheelchair
- Heavy-duty welded construction
- W8 X L7 X H148cm / Prong 19cm

REF CODE	DESCRIPTION / SIZE	UOM	QTY/PACK	MDR
9-IC-IVP201	IV Stand for Wheelchair, 2 Prongs	Ea	1	—





Aluminum Ramp for Wheelchair

Can carry a maximum
load of 272kgs.




Easy to Mount



Anti Skid Surface



Available in 5ft and 6ft.



Easy to Carry



So hurry, book your orders now!



www.inmed.com.ph



Pay by cash, cheque or hassle-free payment options:

Online and over-the-counter payment through banks and credit cards.

INMED CORPORATION

5 Calle Industria, Bagumbayan, Quezon City 1110, Philippines

☎ 02.85711888 📞 📞 **Globe** 0917.8146633

Follow us

